

Product Datasheet

Fas Ligand/TNFSF6 Antibody (14C2) [DyLight 755] NB100-63072IR-0.1ml

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB100-63072IR

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB100-63072IR



NB100-63072IR-0.1ml

Fas Ligand/TNFSF6 Antibody (14C2) [DyLight 755]

| Product Information | |
|----------------------------|---|
| Unit Size | 0.1 ml |
| Concentration | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage | Store at 4C in the dark. |
| Clonality | Monoclonal |
| Clone | 14C2 |
| Preservative | 0.05% Sodium Azide |
| Isotype | IgG1 |
| Conjugate | DyLight 755 |
| Purity | Protein G purified |
| Buffer | 50mM Sodium Borate |

| Product Description | |
|--------------------------------|--|
| Description | This conjugate is made on demand. Actual recovery may vary from the stated volume of this product. The volume will be greater than or equal to the unit size stated on the datasheet. |
| Host | Mouse |
| Gene ID | 356 |
| Gene Symbol | FASLG |
| Species | Human |
| Specificity/Sensitivity | NB100-63072 recognizes the human CD178 cell surface antigen, a 40kDa glycoprotein also known as Fas ligand (CD95L). CD178 is a member of the TNF family, which is expressed by activated T lymphocytes and NK cells. The protein may exist as either a membrane bound or a cleaved soluble form. CD178 plays an important role in T cell development cytotoxicity. Binding of CD178 to Fas (CD95) results in the induction of apoptosis. Clone 14C2 is not reported to have a blocking function. |
| Immunogen | Made to Human CD178 |
| Notes | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries. |

| Product Application Details | |
|------------------------------------|---|
| Applications | ELISA, Flow Cytometry, Immunoprecipitation, Western Blot (Negative) |
| Recommended Dilutions | Flow Cytometry, ELISA, Immunoprecipitation, Western Blot (Negative) |
| Application Notes | Results maybe enhanced using membrane permeabilisation in Flow Cytometry. |





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB100-63072IR-0.1ml

| | |
|---------------|--|
| NBP1-97005IR | Mouse IgG1 Isotype Control (MG1) [DyLight 755] |
| NBP2-49160PEP | Fas Ligand/TNFSF6 Recombinant Protein Antigen |
| 210-TA-005 | TNF-alpha [Unconjugated] |
| DFL00B | Fas Ligand/TNFSF6 [HRP] |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB100-63072IR

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

