

Product Datasheet

Integrin alpha 6/CD49f Antibody (450-30A) [FITC] NB100-63902

Unit Size: 0.1 mg

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB100-63902

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB100-63902



NB100-63902


Integrin alpha 6/CD49f Antibody (450-30A) [FITC]

Product Information	
Unit Size	0.1 mg
Concentration	0.1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	450-30A
Preservative	0.09% Sodium Azide
Isotype	IgG1
Conjugate	FITC
Purity	Protein A purified
Buffer	PBS
Product Description	
Host	Mouse
Gene ID	3655
Gene Symbol	ITGA6
Species	Human
Specificity/Sensitivity	Recognizes CD49f, also known as VLA-6 and alpha 6 integrin. CD49f is expressed by platelets, weakly by monocytes and by a subset of lymphocytes. CD49f is also widely expressed on epithelial tissues.
Immunogen	alpha 6 beta 4 integrin purified from A431 cells
Product Application Details	
Applications	Flow Cytometry
Recommended Dilutions	Flow Cytometry 1:10-1:1000




Images


Western Blot: Integrin alpha 6/CD49f Antibody (450-30A) [FITC] [NB100-63902] - Analysis of Jurkat whole cell lysate.




Flow Cytometry: Integrin alpha 6/CD49f Antibody (450-30A) [FITC] [NB100-63902] - Staining of human peripheral blood monocytes.



Flow Cytometry: Integrin alpha 6/CD49f Antibody (450-30A) [FITC] [NB100-63902] - Staining of human peripheral blood lymphocytes with MOUSE ANTI HUMAN CD49f:FITC



Flow Cytometry: Integrin alpha 6/CD49f Antibody (450-30A) [FITC] [NB100-63902] - Staining of human peripheral blood lymphocytes.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB100-63902

NBL1-12066	Integrin alpha 6/CD49f Overexpression Lysate
NBP1-96849	Mouse IgG1 Isotype Control (MG1) [FITC]
H00003655-Q01-10ug	Recombinant Human Integrin alpha 6/CD49f GST (N-Term) Protein
210-TA-005	TNF-alpha [Unconjugated]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB100-63902

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

