

Product Datasheet

Aminopeptidase N/CD13 Antibody (ER-BMDM1) [Allophycocyanin] NB100-64843APC

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB100-64843APC

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB100-64843APC



NB100-64843APC

Aminopeptidase N/CD13 Antibody (ER-BMDM1) [Allophycocyanin]

| Product Information | |
|---------------------|---|
| Unit Size | 0.1 ml |
| Concentration | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage | Store at 4C in the dark. |
| Clonality | Monoclonal |
| Clone | ER-BMDM1 |
| Preservative | 0.05% Sodium Azide |
| Isotype | IgG2a |
| Conjugate | Allophycocyanin |
| Purity | Protein G purified |
| Buffer | PBS |

| Product Description | |
|-------------------------|---|
| Host | Rat |
| Gene ID | 290 |
| Gene Symbol | ANPEP |
| Species | Mouse |
| Specificity/Sensitivity | NB100-64843 recognizes the murine homologue of human CD13. |
| Immunogen | Aminopeptidase N/CD13 Antibody (ER-BMDM1) made from Balb/c bone marrow derived macrophages. |

| Product Application Details | |
|-----------------------------|--|
| Applications | Flow Cytometry, Immunohistochemistry, Immunohistochemistry-Frozen, Immunohistochemistry-Paraffin |
| Recommended Dilutions | Flow Cytometry, Immunohistochemistry, Immunohistochemistry-Paraffin, Immunohistochemistry-Frozen |
| Application Notes | Optimal dilution of this antibody should be experimentally determined. |

Images

Aminopeptidase N/CD13 Antibody (ER-BMDM1) [Allophycocyanin] [NB100-64843APC] - Vial of APC conjugated antibody. APC is optimally excited at 650 nm by the Red laser (633 or 640 nm) and has an emission maximum of 660 nm.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB100-64843APC

| | |
|---------------|---|
| IC006A | Rat IgG2a Isotype Control (54447) [Allophycocyanin] |
| NB100-64843UV | Aminopeptidase N/CD13 Antibody (ER-BMDM1) [DyLight 350] |
| NBP2-33855PEP | Aminopeptidase N/CD13 Recombinant Protein Antigen |
| 210-TA-005 | TNF-alpha [Unconjugated] |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB100-64843APC

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

