

Product Datasheet

CD43/Sialophorin Antibody (W3/13) - BSA Free NB100-64992-250ug

Unit Size: 250 ug

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 3

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB100-64992

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB100-64992



NB100-64992-250ug


CD43/Sialophorin Antibody (W3/13) - BSA Free

Product Information	
Unit Size	250 ug
Concentration	1.0 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	W3/13
Preservative	0.02% Sodium Azide
Isotype	IgG1
Purity	Protein G purified
Buffer	PBS
Product Description	
Host	Mouse
Gene ID	6693
Gene Symbol	SPN
Species	Rat
Reactivity Notes	Reacts with Rat.
Specificity/Sensitivity	Rat CD43 cell surface antigen.
Immunogen	Rat thymocyte membrane glycoproteins
Product Application Details	
Applications	Western Blot, Flow Cytometry, Immunohistochemistry, Immunohistochemistry-Frozen, Immunohistochemistry-Paraffin
Recommended Dilutions	Western Blot 0.5 ug/ml, Flow Cytometry 1:25-1:100, Immunohistochemistry 5 ug/ml, Immunohistochemistry-Paraffin 5 ug/ml, Immunohistochemistry-Frozen 1:10-1:500
Application Notes	This product requires antigen retrieval using heat treatment prior to staining of paraffin sections. Use in Western blot reported in scientific literature (PMID 7504163).




Images


Western Blot: CD43/Sialophorin Antibody (W3/13) [NB100-64992] - Western blot shows lysates of rat lymph node, rat spleen tissue, and rat thymus tissue. PVDF membrane was probed with 0.5 ug/mL of Mouse Anti-Rat CD43 Monoclonal Antibody (Catalog # NB100-64992) followed by HRP-conjugated Anti-Mouse IgG Secondary Antibody (Catalog # HAF018). A specific band was detected for CD43 at approximately 90 kDa (as indicated). This experiment was conducted under reducing conditions.



Immunohistochemistry-Paraffin: CD43/Sialophorin Antibody (W3/13) [NB100-64992] - CD43 was detected in immersion fixed paraffin-embedded sections of rat thymus using Mouse Anti-Rat CD43 Monoclonal Antibody (Catalog # NB100-64992) at 5 ug/mL for 1 hour at room temperature followed by incubation with the Anti-Mouse IgG VisUCyte™ HRP Polymer Antibody (Catalog # VC001). Before incubation with the primary antibody, tissue was incubated with VisUBlock™ Mouse on Mouse Blocking Reagent (Catalog # VB001) and then subjected to heat-induced epitope retrieval using Antigen Retrieval Reagent-Basic (Catalog # CTS013). Tissue was stained using DAB (brown) and counterstained with hematoxylin (blue). Specific staining was localized to lymphocytes.



Flow Cytometry: CD43/Sialophorin Antibody (W3/13) [NB100-64992] - Staining of rat splenocytes using CD43 FITC.



Publications

Crook K, Hunt SV. Enrichment of early fetal-liver hemopoietic stem cells of the rat using monoclonal antibodies against the transferrin receptor, Thy-1, and MRC-OX82. *Dev Immunol.* [PMID: 8924759] (FLOW, Rat)

Nagai Y, Inobe M, Kikuchi K, Uede T. Functional and phenotypical analysis of subsets of rat CD4+ T cells. *Microbiol Immunol.* [PMID: 7504163] (WB, Rat)

Yamanaka Y, Shigetani Y, Yoshiba K et al. Immunohistochemical analysis of subcutaneous tissue reactions to methacrylate resin-based root canal sealers. *Int Endod J.* 2011-07-01 [PMID: 21375543]



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB100-64992-250ug

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-97005-0.5mg	Mouse IgG1 Isotype Control (MG1)
NB100-64992UV	CD43/Sialophorin Antibody (W3/13) [DyLight 350]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB100-64992

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

