

Product Datasheet

Glycophorin A Antibody (YTH89.1) - BSA Free NB100-65015

Unit Size: 0.1 mg

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB100-65015

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB100-65015



NB100-65015

Glycophorin A Antibody (YTH89.1) - BSA Free

Product Information	
Unit Size	0.1 mg
Concentration	1.0 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	YTH89.1
Preservative	0.09% Sodium Azide
Isotype	IgG2b
Purity	Protein G purified
Buffer	PBS
Product Description	
Description	Novus Biologicals Rat Glycophorin A Antibody (YTH89.1) - BSA Free (NB100-65015) is a monoclonal antibody validated for use in IHC and Flow. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rat
Gene ID	2993
Gene Symbol	GYPA
Species	Human
Specificity/Sensitivity	Recognizes glycophorin A, a major sialoglycoprotein of the human erythrocyte membrane.
Immunogen	Made to Human CD235a
Product Application Details	
Applications	Immunohistochemistry-Paraffin, Flow Cytometry, Immunohistochemistry, Immunohistochemistry-Frozen
Recommended Dilutions	Flow Cytometry 1:10-1:25, Immunohistochemistry 1:10-1:500, Immunohistochemistry-Paraffin 1:100-1:1000, Immunohistochemistry-Frozen 1:100-1:1000
Application Notes	This product does not require antigen retrieval using heat treatment prior to staining of paraffin sections.




Images


Immunohistochemistry-Paraffin: Glycophorin A Antibody (YTH89.1) [NB100-65015] - Human hypophysis showing erythrocyte labelling in sinusoids.



Flow Cytometry: Glycophorin A Antibody (YTH89.1) [NB100-65015] - Analysis using the FITC conjugate of NB100-65015. Staining of human peripheral blood erythrocytes.



Flow Cytometry: Glycophorin A Antibody (YTH89.1) [NB100-65015] - Staining of human peripheral blood erythrocytes with RAT ANTI HUMAN CD235a: RPE.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB100-65015

HAF005	Goat anti-Rat IgG Secondary Antibody [HRP]
F0105B	Goat anti-Rat IgG Secondary Antibody [Phycoerythrin]
DDXCR03	Rat IgG2b Isotype Control
NBP2-59660-50ug	Recombinant Human Glycophorin A His Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB100-65015

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

