

Product Datasheet

CD161 Antibody (B199.2) - BSA Free NB100-65298

Unit Size: 0.1 mg

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB100-65298

Updated 2/21/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB100-65298



NB100-65298


CD161 Antibody (B199.2) - BSA Free

Product Information	
Unit Size	0.1 mg
Concentration	1.0 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	B199.2
Preservative	0.09% Sodium Azide
Isotype	IgG1
Purity	Protein A purified
Buffer	PBS
Product Description	
Host	Mouse
Gene ID	3820
Gene Symbol	KLRB1
Species	Human
Specificity/Sensitivity	CD161 (B199.2)
Immunogen	Purified human NK cells cultured in IL-2
Product Application Details	
Applications	Flow Cytometry, Immunohistochemistry, Immunohistochemistry-Frozen, Immunoprecipitation, Immunofluorescence, Western Blot (Negative)
Recommended Dilutions	Flow Cytometry 1:50-1:100, Immunohistochemistry 1:10-1:500, Immunoprecipitation 1:10-1:500, Immunohistochemistry-Frozen 1:10-1:500, Western Blot (Negative), Immunofluorescence




Images

Flow Cytometry: CD161 Antibody (B199.2) [NB100-65298] - Figure A. PE-A750 conjugated mouse anti human CD3 and Pacific Blue conjugated mouse IgG1 isotype control. Figure B. PE-A750 conjugated mouse anti human CD3 and Pacific Blue conjugated mouse anti human CD161. All experiments performed on red cell lysed human blood gated on lymphocytes in the presence of Human Seroblock.



Flow Cytometry: CD161 Antibody (B199.2) [NB100-65298] - Figure A. PE-A750 conjugated mouse anti human CD3 and A647 conjugated mouse IgG1 isotype control. Figure B. PE-750 conjugated mouse anti human CD3 and A647 conjugated mouse anti human CD161. All experiments performed on red cell lysed human blood gated on lymphocytes in the presence of Human Seroblock.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB100-65298

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-97005-0.5mg	Mouse IgG1 Isotype Control (MG1)
NB100-65298UV	CD161 Antibody (B199.2) [DyLight 350]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB100-65298

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

