

Product Datasheet

Canine Parvovirus Antibody (8H7) - BSA Free NB100-73204

Unit Size: 0.2 mg

Store at 4C. Do not freeze.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB100-73204

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB100-73204



NB100-73204**Canine Parvovirus Antibody (8H7) - BSA Free**

| Product Information | |
|------------------------------------|---|
| Unit Size | 0.2 mg |
| Concentration | Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services. |
| Storage | Store at 4C. Do not freeze. |
| Clonality | Monoclonal |
| Clone | 8H7 |
| Preservative | 0.09% Sodium Azide |
| Isotype | IgG2a |
| Purity | Protein A purified |
| Buffer | PBS (pH 7.4) |
| Product Description | |
| Description | Novus Biologicals Mouse Canine Parvovirus Antibody (8H7) - BSA Free (NB100-73204) is a monoclonal antibody validated for use in Immunodiffusion, Haemagglutination, Sandwich ELISA, Lateral Flow Assay. All Novus Biologicals antibodies are covered by our 100% guarantee. |
| Host | Mouse |
| Species | Virus |
| Immunogen | Hybridoma clone has been derived from hybridization of Sp2/0 myeloma cells with spleen cells of Balb/c mice immunized with Canine parvovirus. |
| Product Application Details | |
| Applications | Haemagglutination, Immunodiffusion, Sandwich ELISA, Lateral Flow Assay |
| Recommended Dilutions | Immunodiffusion, Haemagglutination, Sandwich ELISA 1:100-1:2000, Lateral Flow Assay |
| Application Notes | Recommended pair for sandwich immunoassay (capture-detection): NB110-2521 - NB100-73204 |



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB100-73204

| | |
|------------|--|
| HAF007 | Goat anti-Mouse IgG Secondary Antibody [HRP] |
| NB7539 | Goat anti-Mouse IgG (H+L) Secondary Antibody [HRP] |
| NBP1-96778 | Mouse IgG2a Isotype Control (M2A) |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB100-73204

Earn gift cards/discounts by submitting a publication using this product:

www.novusbio.com/publications

