

Product Datasheet

Melanopsin Antibody NB100-74460

Unit Size: 100 ug

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB100-74460

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB100-74460



NB100-74460

Melanopsin Antibody

Product Information

Unit Size	100 ug
Concentration	1 mg/ml
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.05% Sodium Azide
Isotype	IgG
Purity	Immunogen affinity purified
Buffer	PBS, 1 mg/mL BSA

Product Description


Host	Rabbit
Gene ID	94233
Gene Symbol	OPN4
Species	Human, Mouse, Rat
Reactivity Notes	Human reactivity reported in scientific literature (PMID: 20971788). Mouse (PMID: 22886938).
Immunogen	Synthetic peptide corresponding to residues E(455) Q K S K T P K T K R H L P S L D R R M(474) of rat melanopsin.

Product Application Details

Applications	Western Blot, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Frozen, Knockout Validated
Recommended Dilutions	Western Blot 0.1 - 0.2 ug/mL, Immunohistochemistry 1:10 - 1:500, Immunocytochemistry/ Immunofluorescence 5 ug/mL, Immunohistochemistry-Frozen 1:1000, Knockout Validated
Application Notes	Successful WB, IF, and ICC usages. WB: Detects approx. 53 and 85 kDa proteins representing overexpressed unglycosylated and glycosylated melanopsin from HEK 293 cells (does not work well on endogenously expressed protein). ICC/IF: NB100-74460 detects melanopsin from rat retinal ganglion cells (RGCs). Knockout data (PMID: 31924821).

Images

Immunocytochemistry/Immunofluorescence: Melanopsin Antibody [NB100-74460] - Staining in rat retinal ganglion cells.

Fig. 1

Publications

Roy D, Levi K, Kiss V et al. Rhodopsin and melanopsin coexist in mammalian sperm cells and activate different signaling pathways for thermotaxis Sci Rep 2020-01-10 [PMID: 31924821] (KO, ICC/IF, Mouse, Human)





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB100-74460

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
664-LI-025	LIGHT/TNFSF14 [Unconjugated]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB100-74460

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

