

Product Datasheet

TIMP-2 Antibody (F27 P3 A4)

NB100-74552

Unit Size: 100ug

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB100-74552

Updated 4/29/2021 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB100-74552



NB100-74552

TIMP-2 Antibody (F27 P3 A4)


Product Information	
Unit Size	100ug
Concentration	1 mg/ml
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	F27 P3 A4
Preservative	0.05% Sodium Azide
Isotype	IgG2a Kappa
Purity	Protein A purified
Buffer	PBS with 1 mg/ml BSA

Product Description	
Host	Mouse
Gene ID	7077
Gene Symbol	TIMP2
Species	Human, Mouse
Reactivity Notes	Please note that this antibody is reactive to Mouse and derived from the same host, Mouse. Additional Mouse on Mouse blocking steps may be required for IHC and ICC experiments. Please contact Technical Support for more information.
Immunogen	Synthetic peptide corresponding to residues D(56) S G N D I Y G N P I(66) of mouse TIMP-2.


Product Application Details	
Applications	Western Blot, Immunocytochemistry/Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin
Recommended Dilutions	Western Blot 1:500, Immunohistochemistry 1:100 - 1:500, Immunocytochemistry/Immunofluorescence 1 - 2 ug/ml, Immunohistochemistry-Paraffin 1:50 - 1:500
Application Notes	WB: Detects an approx. 18-24 kDa protein representing TIMP-2 from mouse platelets. It is recommended to block with 3% BSA.

Images


Immunocytochemistry/Immunofluorescence: TIMP-2 Antibody (F27 P3 A4) [NB100-74552] - Cells were grown on chamber slides and fixed with formaldehyde prior to staining. Cells were probed without (control) or with a TIMP2 monoclonal antibody at a dilution of 1:200 overnight at 4 C, washed with PBS and incubated with a DyLight-488 conjugated secondary antibody and nuclei with DAPI (blue) is shown.




Immunohistochemistry-Paraffin: TIMP-2 Antibody (F27 P3 A4) [NB100-74552] - Normal biopsies of deparaffinized Mouse pancreas tissue at a dilution of 1:100.




Immunocytochemistry/Immunofluorescence: TIMP-2 Antibody (F27 P3 A4) [NB100-74552] - Cells were grown on chamber slides and fixed with formaldehyde prior to staining. Cells were probed without (control) or with a TIMP2 monoclonal antibody at a dilution of 1:200 overnight at 4 C, washed with PBS and incubated with a DyLight-488 conjugated secondary antibody and nuclei with DAPI (blue) is shown.



Immunocytochemistry/Immunofluorescence: TIMP-2 Antibody (F27 P3 A4) [NB100-74552] - Cells were grown on chamber slides and fixed with formaldehyde prior to staining. Cells were probed without (control) or with a TIMP2 monoclonal antibody at a dilution of 1:200 overnight at 4 C, washed with PBS and incubated with a DyLight-488 conjugated secondary antibody and nuclei with DAPI (blue) is shown.



Immunohistochemistry-Paraffin: TIMP-2 Antibody (F27 P3 A4) [NB100-74552] - Normal biopsies of deparaffinized Human pancreas tissue at a dilution of 1:500.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA

Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
novus@novusbio.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada

Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom

Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: technical@novusbio.com
Orders: orders@novusbio.com
General: novus@novusbio.com

Products Related to NB100-74552

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP (Horseradish Peroxidase)]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-96981-0.5mg	Mouse IgG2a Kappa Isotype Control (M2AK)
NBP2-61356-10ug	Recombinant Human TIMP-2 Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB100-74552

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

