

Product Datasheet

MYO18A Antibody - BSA Free

NB100-79839

Unit Size: 100 ul

Store at 4C. Do not freeze.

www.novusbio.com



technical@novusbio.com

Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB100-79839

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB100-79839



NB100-79839

MYO18A Antibody - BSA Free

Product Information

Unit Size	100 ul
Concentration	1.0 mg/ml
Storage	Store at 4C. Do not freeze.
Clonality	Polyclonal
Preservative	0.09% Sodium Azide
Isotype	IgG
Purity	Immunogen affinity purified
Buffer	Tris-Citrate/Phosphate (pH 7.0 - 8.0)

Product Description


Description	Novus Biologicals Rabbit MYO18A Antibody - BSA Free (NB100-79839) is a polyclonal antibody validated for use in IP. Anti-MYO18A Antibody: Cited in 1 publication. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	399687
Gene Symbol	MYO18A
Species	Human
Immunogen	The immunogen this antibody was made to, maps to a region between residue 300 and 350 of human myosin 18A using the numbering given in entry NP_510880.2 (GeneID 399687).

Product Application Details

Applications	Immunoprecipitation
Recommended Dilutions	Immunoprecipitation 2-5 mcg/mg lysate
Application Notes	This antibody has not performed satisfactorily when used for WB of Myo18A/Myosin XVIII A in crude preparations (e.g. whole cell lysate). This antibody can be used for WB of enriched (e.g. immunoprecipitated) sources of Myo18A/Myosin XVIII A.

Images

Immunoprecipitation: MYO18A Antibody [NB100-79839] - Detection of Human Myo18A/Myosin XVIII A by Western Blot of Immunoprecipitates. Samples: RIPA whole cell lysate (1 mg for IP, 20% of IP loaded) from HeLa cells. Antibodies: Affinity purified rabbit anti-Myo18A/Myosin XVIII A antibody NB100-79839 used for IP at 3 mcg/mg lysate. Myo18A/Myosin XVIII A was also immunoprecipitated by rabbit anti- Myo18A/Myosin XVIII A antibody NB100-79840, which recognizes a downstream epitope. Detection: Chemiluminescence with an exposure time of 1 second.



Publications

Aguilar-Cuenca R, Llorente-Gonzalez C, Chapman JR et al. Tyrosine Phosphorylation of the Myosin Regulatory Light Chain Controls Non-muscle Myosin II Assembly and Function in Migrating Cells Curr. Biol. 2020-06-02 [PMID: 32502416]





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB100-79839

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NBP1-83021PEP	MYO18A Recombinant Protein Antigen

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB100-79839

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

