

Product Datasheet

IRS1 [p Ser312] Antibody - BSA Free NB100-81999

Unit Size: 0.05 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB100-81999

Updated 2/21/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB100-81999



NB100-81999


IRS1 [p Ser312] Antibody - BSA Free

| Product Information | |
|------------------------------------|---|
| Unit Size | 0.05 ml |
| Concentration | 1.0 mg/ml |
| Storage | Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles. |
| Clonality | Polyclonal |
| Preservative | 0.02% Sodium Azide |
| Isotype | IgG |
| Purity | Immunogen affinity purified |
| Buffer | PBS (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150 mM NaCl, 50% glycerol |
| Target Molecular Weight | 170 kDa |
| Product Description | |
| Host | Rabbit |
| Gene ID | 3667 |
| Gene Symbol | IRS1 |
| Species | Human, Mouse, Rat |
| Specificity/Sensitivity | Detects endogenous levels of IRS1 only when phosphorylated at serine 312. |
| Immunogen | The antiserum was produced against synthesized phosphopeptide derived from human IRS1 around the phosphorylation site of serine 312 (A-T-SP-P-A). |
| Product Application Details | |
| Applications | Western Blot |
| Recommended Dilutions | Western Blot 1:500-1:1000 |




Images

Western Blot: IRS1 [p Ser312] Antibody [NB100-81999] - Analysis of extract from 293 cells untreated or treated with insulin using IRS-1 (Phospho-Ser312) antibody



Western Blot: IRS1 [p Ser312] Antibody [NB100-81999] - Western blot analysis of extract from 293 cells, using IRS-1 (Phospho-Ser312) antibody (Lane 1 and 2).





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB100-81999

| | |
|---------------|---|
| HAF008 | Goat anti-Rabbit IgG Secondary Antibody [HRP] |
| NB7160 | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| NBP2-24891 | Rabbit IgG Isotype Control |
| NBP2-68666PEP | IRS1 Recombinant Protein Antigen |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB100-81999

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

