

Product Datasheet

PRPF4B Antibody NB100-86997

Unit Size: 0.1 ml

Store at 4C. Do not freeze.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB100-86997

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB100-86997



NB100-86997**PRPF4B Antibody**


| Product Information | |
|----------------------------|-----------------------------|
| Unit Size | 0.1 ml |
| Concentration | 0.2 mg/ml |
| Storage | Store at 4C. Do not freeze. |
| Clonality | Polyclonal |
| Preservative | 0.09% Sodium Azide |
| Isotype | IgG |
| Purity | Immunogen affinity purified |
| Buffer | TBS and 0.1% BSA |

| Product Description | |
|----------------------------|--|
| Host | Rabbit |
| Gene ID | 8899 |
| Gene Symbol | PRPF4B |
| Species | Human |
| Reactivity Notes | Orangutan (100%). |
| Immunogen | The immunogen recognized by this antibody maps to a region between residue 1 and 50 of human PRP4 pre-mRNA processing factor 4 homolog B using the numbering given in entry NP_003904.3 (GenelD 8899). |


| Product Application Details | |
|------------------------------------|--|
| Applications | Western Blot, Immunohistochemistry, Immunohistochemistry-Paraffin, Immunoprecipitation (Negative) |
| Recommended Dilutions | Western Blot 1:2000-1:10000, Immunohistochemistry 1:200 -1:1000, Immunohistochemistry-Paraffin 1:200 -1:1000, Immunoprecipitation (Negative) |
| Application Notes | Epitope retrieval with citrate buffer pH6.0 is recommended for FFPE tissue sections. Immunoprecipitation Not recommended. |

Images

Western Blot: PRPF4B Antibody [NB100-86997] - Whole cell lysate from HeLa (5, 15 and 50 ug) and 293T (T; 50 ug) cells. Antibody used at 0.04 ug/ml.



Immunohistochemistry-Paraffin: PRPF4B Antibody [NB100-86997] - Human breast carcinoma. Antibody: Affinity purified rabbit anti- PRPF4B used at a dilution of 1:200 (1ug/ml). Detection: DAB





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB100-86997

| | |
|------------|---|
| NBL1-14823 | PRPF4B Overexpression Lysate |
| HAF008 | Goat anti-Rabbit IgG Secondary Antibody [HRP] |
| NB7160 | Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP] |
| NBP2-24891 | Rabbit IgG Isotype Control |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB100-86997

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

