

Product Datasheet

VACHT/SLC18A3 Antibody NB100-91348

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 2

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB100-91348

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB100-91348



NB100-91348**VACHT/SLC18A3 Antibody**


Product Information	
Unit Size	0.1 ml
Concentration	This product is unpurified. The exact concentration of antibody is not quantifiable.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	No Preservative
Reconstitution Instructions	Reconstitute in 0.1 ml of sterile water. Centrifuge to remove any insoluble material. Glycerol may be added (1:1) for additional stability. Please note the sample size is provided in reconstituted format.
Isotype	IgG
Purity	Unpurified
Buffer	Lyophilized from whole antisera

Product Description	
Host	Rabbit
Gene ID	6572
Gene Symbol	SLC18A3
Species	Human, Mouse, Rat, Monkey
Reactivity Notes	Marmoset
Immunogen	A synthetic peptide from human VACHT/SLC18A3 conjugated to blue carrier protein was used as the antigen. The antigen is homologous in rat and mouse.


Product Application Details	
Applications	Western Blot, Immunohistochemistry, Immunohistochemistry-Paraffin
Recommended Dilutions	Western Blot 1:2000, Immunohistochemistry 1:1000, Immunohistochemistry-Paraffin 1:1000

Images

Western Blot: VACHT/SLC18A3 Antibody [NB100-91348] - Rat and mouse tissue lysates. Blocking: 1% LFDM for 30 min at RT; primary antibody: dilution 1:2000 incubated at 4C overnight.




Immunohistochemistry-Paraffin: VACHT/SLC18A3 Antibody [NB100-91348] - Rat spinal cord. The animal was perfused at a pressure of 130 mmHg with 300 ml 4% FA before being processed for paraffin embedding. HIER: Tris-EDTA, pH 9 for 20 min. Blocking: 0.2% LFDM in TBST filtered thru 0.2 μ m.




Immunohistochemistry-Paraffin: VACHT/SLC18A3 Antibody [NB100-91348] - IHC on rat spinal cord (paraffin section) using Rabbit antibody to (Vesicular Acetylcholine Transporter) VACHT: whole serum at 1 : 200 dilution incubated overnight at 4C, developed with DAB/Ni and counterstained with hematoxylin.




Immunohistochemistry-Paraffin: VACHT/SLC18A3 Antibody [NB100-91348] - IHC on rat spinal cord (paraffin section) using Rabbit antibody to (Vesicular Acetylcholine Transporter) VACHT: whole serum at 1 : 200 dilution incubated overnight at 4C, developed with DAB/Ni and counterstained with hematoxylin.




Immunohistochemistry-Paraffin: VACHT/SLC18A3 Antibody [NB100-91348] - IHC on rat spinal cord (paraffin section) using Rabbit antibody to (Vesicular Acetylcholine Transporter) VACHT: whole serum at 1 : 200 dilution incubated overnight at 4C, developed with DAB/Ni and counterstained with hematoxylin.




Immunohistochemistry-Paraffin: VACHT/SLC18A3 Antibody [NB100-91348] - IHC on rat spinal cord (paraffin section) using Rabbit antibody to (Vesicular Acetylcholine Transporter) VACHT: whole serum at 1 : 200 dilution incubated overnight at 4C, developed with DAB/Ni and counterstained with hematoxylin.



Immunohistochemistry-Paraffin: VACHT/SLC18A3 Antibody [NB100-91348] - IHC on rat spinal cord (paraffin section) using Rabbit antibody to (Vesicular Acetylcholine Transporter) VACHT: whole serum at 1 : 200 dilution incubated overnight at 4C, developed with DAB/Ni and counterstained with hematoxylin.



Immunohistochemistry-Paraffin: VACHT/SLC18A3 Antibody [NB100-91348] - IHC on rat spinal cord (paraffin section) using Rabbit antibody to (Vesicular Acetylcholine Transporter) VACHT: whole serum at 1 : 200 dilution incubated overnight at 4C, developed with DAB/Ni and counterstained with hematoxylin.



Publications

Tie P, Cheng J, Xue MX et al. SLC18A3 promoted renal cancer development through acetylcholine/cAMP signaling American journal of cancer research 2022-09-15 [PMID: 36225635] (IHC-P, WB, Human)

Gokcek-Sarac C, Akcay G, Karakurt S et al. Possible effects of different doses of 2.1 GHz electromagnetic radiation on learning, and hippocampal levels of cholinergic biomarkers in Wistar rats Electromagnetic biology and medicine 2020-12-01 [PMID: 33259237] (WB)



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB100-91348

NB100-91348PEP	VACHT/SLC18A3 Antibody Blocking Peptide
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB100-91348

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

