

Product Datasheet

RanBP2 Antibody NB100-93337

Unit Size: 0.1 ml

Store at 4C. Do not freeze.

www.novusbio.com



technical@novusbio.com

Publications: 1


Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB100-93337

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB100-93337



NB100-93337

RanBP2 Antibody

Product Information	
Unit Size	0.1 ml
Concentration	0.2 mg/ml
Storage	Store at 4C. Do not freeze.
Clonality	Polyclonal
Preservative	0.09% Sodium Azide
Isotype	IgG
Purity	Immunogen affinity purified
Buffer	TBS and 0.1% BSA
Target Molecular Weight	358 kDa


Product Description	
Description	Novus Biologicals Rabbit RanBP2 Antibody (NB100-93337) is a polyclonal antibody validated for use in WB and IP. Anti-RanBP2 Antibody: Cited in 1 publication. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	5903
Gene Symbol	RANBP2
Species	Human, Mouse
Reactivity Notes	Chimpanzee (100%).
Immunogen	The immunogen recognized by this antibody maps to a region between residue 1050 and 1100 of human RAN binding protein 2 using the numbering given in entry NP_006258.3 (GenelD 5903).

Product Application Details	
Applications	Western Blot, Immunoprecipitation
Recommended Dilutions	Western Blot 1:2000-1:10000, Immunoprecipitation 2 - 10 ug/mg lysate




Images


Western Blot: RanBP2 Antibody [NB100-93337] - Whole cell lysate (50 ug) from NIH3T3 and TCMK-1 cells prepared using NETN lysis buffer. Antibody: Affinity purified rabbit anti-RanBP2 antibody used for WB at 1 ug/ml. Detection: Chemiluminescence with an exposure time of 30 seconds.



Western Blot: RanBP2 Antibody [NB100-93337] - Detection of Human RanBP2 on HeLa whole cell lysate using NB100-93337. RanBP2 was also immunoprecipitated by rabbit anti-RanBP2 antibody NB100-93338.



Western Blot: RanBP2 Antibody [NB100-93337] - Detection of Human RanBP2 by Western Blot. Samples: Whole cell lysate (50 ug) from HeLa and 293T cells prepared using NETN lysis buffer. Antibody: Affinity purified rabbit anti-RanBP2 antibody NB100-93337 used for WB at 0.1 ug/ml. Detection: Chemiluminescence with an exposure time of 3 minutes.



Publications

Wang WF, Yan L, Liu Z et al. HSP70-Hrd1 axis precludes the oncorepressor potential of N-terminal misfolded Blimp-1s in lymphoma cells Nat Commun 2017-08-25 [PMID: 28842558] (WB, Mouse)



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB100-93337

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB100-93337

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

