

Product Datasheet

ADNP Antibody NB110-40556

Unit Size: 0.1 ml

Store at 4C. Do not freeze.

www.novusbio.com



technical@novusbio.com

Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB110-40556

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB110-40556



NB110-40556**ADNP Antibody**

Product Information	
Unit Size	0.1 ml
Concentration	0.25 mg/ml
Storage	Store at 4C. Do not freeze.
Clonality	Polyclonal
Preservative	0.09% Sodium Azide
Isotype	IgG
Purity	Immunogen affinity purified
Buffer	TBS and 0.1% BSA


Product Description	
Description	Novus Biologicals Rabbit ADNP Antibody (NB110-40556) is a polyclonal antibody validated for use in IHC. Anti-ADNP Antibody: Cited in 1 publication. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	23394
Gene Symbol	ADNP
Species	Human, Rat
Reactivity Notes	Based on 100% sequence identity, this antibody is predicted to react with Gorilla, White-tufted-ear Marmoset, Crab-eating Macaque, Small-eared galago and Northern White-cheeked Gibbon. Rat reactivity reported in scientific literature (PMID: 18571851)
Immunogen	The immunogen recognized by this antibody maps to a region between residues 1050 and the C-terminus (residue 1102) of human Activity-Dependent Neuroprotective Protein using the numbering given in entry NP_056154.1 (GeneID 23394)

Product Application Details	
Applications	Immunohistochemistry-Paraffin, Immunohistochemistry
Recommended Dilutions	Immunohistochemistry 1:100-1:500, Immunohistochemistry-Paraffin 1:100-1:500




Images

Immunohistochemistry: ADNP Antibody [NB110-40556] - Sample: FFPE section of human metastatic bone marrow. Antibody: Affinity purified rabbit anti-ADNP used at a dilution of 1:250. Detection: DAB staining using Immunohistochemistry Accessory Kit .



Immunohistochemistry: ADNP Antibody [NB110-40556] - Sample: FFPE section of human breast adenocarcinoma. Antibody: Affinity purified rabbit anti-ADNP used at a dilution of 1:250. Detection: DAB staining using Immunohistochemistry Accessory Kit ..



Publications

Cattano D, Valleggi S, Ma D et al. Xenon induces transcription of ADNP in neonatal rat brain. Neurosci Lett 2008-08-08 [PMID: 18571851] (Rat)



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB110-40556

NBL1-07346	ADNP Overexpression Lysate
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB110-40556

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

