

# Product Datasheet

## OXPAT Antibody [PerCP] NB110-60511PCP

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NB110-60511PCP](http://www.novusbio.com/NB110-60511PCP)

Updated 10/31/2023 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NB110-60511PCP](http://www.novusbio.com/reviews/destination/NB110-60511PCP)



**NB110-60511PCP**

OXPAT Antibody [PerCP]

| Product Information         |                                                                                                                        |
|-----------------------------|------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                                 |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                            |
| Storage                     | Store at 4C in the dark.                                                                                               |
| Clonality                   | Polyclonal                                                                                                             |
| Preservative                | 0.05% Sodium Azide                                                                                                     |
| Isotype                     | IgG                                                                                                                    |
| Conjugate                   | PerCP                                                                                                                  |
| Purity                      | Immunogen affinity purified                                                                                            |
| Buffer                      | PBS                                                                                                                    |
| Product Description         |                                                                                                                        |
| Host                        | Rabbit                                                                                                                 |
| Gene ID                     | 440503                                                                                                                 |
| Gene Symbol                 | PLIN5                                                                                                                  |
| Species                     | Human, Mouse                                                                                                           |
| Reactivity Notes            | Immunogen displays the following percentage of sequence identity for non-tested species: Rat (89%).                    |
| Immunogen                   | A synthetic peptide made to a C-terminal region of the mouse OXPAT protein, within residues 400-463. [UniProt# Q8BVZ1] |
| Product Application Details |                                                                                                                        |
| Applications                | Western Blot, Simple Western                                                                                           |
| Recommended Dilutions       | Western Blot, Simple Western                                                                                           |
| Application Notes           | Optimal dilution of this antibody should be experimentally determined.                                                 |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA

Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
novus@novusbio.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada

Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: technical@novusbio.com  
Orders: orders@novusbio.com  
General: novus@novusbio.com

### **Products Related to NB110-60511PCP**

---

|                |                                    |
|----------------|------------------------------------|
| NBP2-24891PCP  | Rabbit IgG Isotype Control [PerCP] |
| NB110-60511PEP | OXPAT Antibody Blocking Peptide    |
| DRP300         | Adiponectin/Acrp30 [HRP]           |
| NB110-60509PEP | OXPAT Antibody Blocking Peptide    |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NB110-60511PCP](http://www.novusbio.com/reviews/submit/NB110-60511PCP)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

