

Product Datasheet

Tenascin C Antibody (4C8MS) [Allophycocyanin/Cy7] NB110-68136APCCY7

Unit Size: 0.1 ml

Store at 4C in the dark. Do not freeze.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB110-68136APCCY7

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB110-68136APCCY7



NB110-68136APCCY7

Tenascin C Antibody (4C8MS) [Allophycocyanin/Cy7]

| Product Information | |
|------------------------------------|---|
| Unit Size | 0.1 ml |
| Concentration | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage | Store at 4C in the dark. Do not freeze. |
| Clonality | Monoclonal |
| Clone | 4C8MS |
| Preservative | 0.05% Sodium Azide |
| Isotype | IgG1 |
| Conjugate | Allophycocyanin/Cy7 |
| Purity | Protein A purified |
| Buffer | PBS |
| Product Description | |
| Host | Mouse |
| Gene ID | 3371 |
| Gene Symbol | TNC |
| Species | Human, Mouse, Rat, Feline |
| Reactivity Notes | Feline reactivity reported in customer review. |
| Specificity/Sensitivity | NB110-68136 specifically reacts with Domain B on FNIII repeats of Tenascin C. |
| Immunogen | Recombinant human Tenascin C [Swiss-Prot# P24821]. |
| Product Application Details | |
| Applications | Western Blot, ELISA, Flow Cytometry, Flow (Intracellular), Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Frozen, Immunohistochemistry-Paraffin, CyTOF-ready |
| Recommended Dilutions | Western Blot, Flow Cytometry, ELISA, Immunohistochemistry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry-Paraffin, Immunohistochemistry-Frozen, Flow (Intracellular), CyTOF-ready |
| Application Notes | Optimal dilution of this antibody should be experimentally determined. |





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB110-68136APCCY7

| | |
|------------------|--|
| NBP1-97005APCCY7 | Mouse IgG1 Isotype Control (MG1) [Allophycocyanin/Cy7] |
| NBP1-89639PEP | Tenascin C Recombinant Protein Antigen |
| 236-EG-200 | EGF [Unconjugated] |
| 3358-TC-050 | Tenascin C [Unconjugated] |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB110-68136APCCY7

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

