

Product Datasheet

Fibroblast Activation Protein alpha/FAP Antibody [PE] NB110-85534PE

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB110-85534PE

Updated 10/6/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB110-85534PE



NB110-85534PE**Fibroblast Activation Protein alpha/FAP Antibody [PE]**


Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark.
Clonality	Polyclonal
Preservative	0.05% Sodium Azide
Isotype	IgG
Conjugate	PE
Purity	Immunogen affinity purified
Buffer	PBS

Product Description	
Host	Rabbit
Gene ID	2191
Gene Symbol	FAP
Species	Human, Mouse
Reactivity Notes	A BLAST analysis was used to suggest partial cross-reactivity with Fibroblast Activation Protein alpha/FAP from human and rat sources based on ~88% homology with the immunizing sequence. Reactivity with Fibroblast Activation Protein alpha/FAP from other sources has not been determined. Human reactivity reported from a verified customer review.
Immunogen	This affinity purified Fibroblast Activation Protein alpha/FAP Antibody was prepared from whole rabbit serum produced by repeated immunizations with a synthetic peptide corresponding to an internal region of mouse Fibroblast Activation Protein alpha/Fibroblast Activation Protein alpha/FAP protein. (Uniprot: P97321)

Product Application Details	
Applications	Western Blot, ELISA, Flow Cytometry
Recommended Dilutions	Western Blot, Flow Cytometry Validated for Flow Cytometry from a verified customer review., ELISA
Application Notes	Optimal dilution of this antibody should be experimentally determined.

Images

Flow Cytometry: Fibroblast Activation Protein alpha/FAP Antibody [PE] [NB110-85534PE] - Specific co-culture system created to study human macrophage cancer cells' conduction of MMT process, and Fibroblast Activation Protein alpha/FAP Antibody [Alexa Fluor® 488] (Catalog # NB110-85534AF488) was used as a myo-CAFs marker. Image from a verified customer review.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA

Phone: 303.730.1950

Toll Free: 1.888.506.6887

Fax: 303.730.1966

nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada

Phone: 905.827.6400

Toll Free: 855.668.8722

Fax: 905.827.6402

canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom

Phone: (44) (0) 1235 529449

Free Phone: 0800 37 34 15

Fax: (44) (0) 1235 533420

info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com

Technical Support: nb-technical@bio-techne.com

Orders: nb-customerservice@bio-techne.com

General: novus@novusbio.com

Products Related to NB110-85534PE

NBP2-24983	Rabbit IgG Isotype Control [PE]
DPC900	Fibroblast Activation Protein alpha/FAP [HRP]
AF1180	DPPIV/CD26 Antibody [Unconjugated]
AF2335	Aminopeptidase N/CD13 Antibody [Unconjugated]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB110-85534PE

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

