

Product Datasheet

Cytokeratin 1 Antibody (3C10F7) - BSA Free NB110-89450

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 2

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB110-89450

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB110-89450




NB110-89450**Cytokeratin 1 Antibody (3C10F7) - BSA Free**


Product Information	
Unit Size	0.1 ml
Concentration	This product is unpurified. The exact concentration of antibody is not quantifiable.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	3C10F7
Preservative	0.03% Sodium Azide
Isotype	IgG1
Purity	Unpurified
Buffer	Ascites
Product Description	
Host	Mouse
Gene ID	3848
Gene Symbol	KRT1
Species	Human, Porcine, Bovine, Canine, Equine, Feline, Goat
Reactivity Notes	Porcine, Bovine, Goat, Equine, Canine, and Feline reactivities reported in scientific literature (PMID: 26190143).
Specificity/Sensitivity	CK1 (PMC01)
Immunogen	Purified recombinant fragment of Cytokeratin 1 expressed in E. Coli.
Product Application Details	
Applications	Western Blot, ELISA, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin, Immunohistochemistry Free-Floating
Recommended Dilutions	Western Blot 1:500 - 1:2000, ELISA 1:10000, Immunohistochemistry 1:10 - 1:500, Immunocytochemistry/ Immunofluorescence 2-5 ug/ml, Immunohistochemistry-Paraffin 1:200-1:1000, Immunohistochemistry Free-Floating
Application Notes	Use in Immunohistochemistry free floating reported in scientific literature (PMID: 26190143).

Images


Western Blot: Cytokeratin 1 Antibody (3C10F7) [NB110-89450] - Analysis using CK1 mouse mAb against truncated CK1 recombinant protein.



Immunocytochemistry/Immunofluorescence: Cytokeratin 1 Antibody (3C10F7) [NB110-89450] - HeLa cells were fixed for 10 minutes using 10% formalin and then permeabilized for 5 minutes using 1X PBS + 0.05% Triton-X100. The cells were incubated with anti-Cytokeratin 1 (3C10F7) at 5 ug/ml overnight at 4C and detected with an anti-mouse Dylight 488 (Green) at a 1:500 dilution. Actin was detected with Phalloidin 568 (Red) at a 1:200 dilution. Nuclei were counterstained with DAPI (Blue). Cells were imaged using a 40X objective.



Immunohistochemistry-Paraffin: Cytokeratin 1 Antibody (3C10F7) [NB110-89450] - Analysis of human liver tissue (A), colon carcinoma (B), lung carcinoma (C) and esophagus tissue (D), showing membrane localization using CK1 mouse mAb with DAB staining.



Publications

Fuchs C, Stalnaker KJ, Dalgard CL et al. Plantar Skin Exhibits Altered Physiology, Constitutive Activation of Wound-Associated Phenotypes, and Inherently Delayed Healing. *The Journal of investigative dermatology* 2024-06-22 [PMID: 38237729]

Nautscher N, Bauer A, Steffl M, Amselgruber WM et al. Comparative morphological evaluation of domestic animal cornea *Vet Ophthalmol.* 2015-07-20 [PMID: 26190143] (IHC-P, Equine, Feline, Goat, Bovine, Canine, Porcine)

Details:

Cytokeratin 1 antibody used for Immunohistochemistry on paraffin embedded corneal sections from pig, cow, goat, horse, dog and cat.



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB110-89450

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-97005-0.5mg	Mouse IgG1 Isotype Control (MG1)
H00003848-Q01-10ug	Recombinant Human Cytokeratin 1 GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB110-89450

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

