

Product Datasheet

Troponin I Type 3 (cardiac) Antibody (MF4cc) [Allophycocyanin] NB120-10239APC

Unit Size: 0.2 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB120-10239APC

Updated 10/30/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB120-10239APC



NB120-10239APC

Troponin I Type 3 (cardiac) Antibody (MF4cc) [Allophycocyanin]

Product Information	
Unit Size	0.2 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark.
Clonality	Monoclonal
Clone	MF4cc
Preservative	0.05% Sodium Azide
Isotype	IgG1
Conjugate	Allophycocyanin
Purity	Protein A purified
Buffer	PBS
Product Description	
Host	Mouse
Gene ID	7137
Gene Symbol	TNNI3
Species	Human, Rat, Porcine, Bovine, Canine, Feline, Goat, Fish (Negative), Rabbit (Negative)
Reactivity Notes	Cross-reacts with Rat, Cow, Cat, Dog, Goat, Pig and Human. Does not cross-react with Fish or Rabbit. Not yet tested in other species.
Specificity/Sensitivity	No cross-reaction with skeletal TnT
Immunogen	Hybridoma clone has been derived from hybridization of Sp2/0 myeloma cells with spleen cells of Balb/c mice immunized with free human cTnI.
Notes	This conjugate is made on demand. Actual recovery may vary from the stated volume of this product. The volume will be greater than or equal to the unit size stated on the datasheet.
Product Application Details	
Applications	Western Blot, ELISA, Immunoassay, Immunohistochemistry, Immunoprecipitation, Sandwich ELISA, Immunoaffinity Purification
Recommended Dilutions	Western Blot, ELISA, Immunohistochemistry, Immunoprecipitation, Immunoassay, Sandwich ELISA, Immunoaffinity Purification
Application Notes	Optimal dilution of this antibody should be experimentally determined.



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB120-10239APC

NBP1-97005APC	Mouse IgG1 Isotype Control (MG1) [Allorhycocyanin]
H00007137-P01-2ug	Recombinant Human Troponin I Type 3 (cardiac) GST (N-Term) Protein
210-TA-005	TNF-alpha [Unconjugated]
NBP3-00456	Mouse Troponin I Type 3 (cardiac) Sandwich ELISA Optimization Kit

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB120-10239APC

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

