

Product Datasheet

NUP62 Antibody (RL31) NB120-2736

Unit Size: 100 ug

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB120-2736

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB120-2736



NB120-2736

NUP62 Antibody (RL31)


Product Information	
Unit Size	100 ug
Concentration	0.5 mg/ml
Storage	Store at -20C. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	RL31
Preservative	0.05% Sodium Azide
Isotype	IgG
Purity	Protein A purified
Buffer	PBS (pH 7.2-7.4) with 1mg/ml BSA and 30% glycerol

Product Description	
Description	Novus Biologicals Armenian Hamster NUP62 Antibody (RL31) (NB120-2736) is a monoclonal antibody validated for use in WB. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Armenian Hamster
Gene ID	23636
Gene Symbol	NUP62
Species	Rat
Specificity/Sensitivity	Detects the nucleoporin p62 O-linked glycoprotein in the p62 complex from rat tissues.
Immunogen	O-linked glycoprotein-enriched fraction from rat liver nuclear envelopes raised in Armenian hamsters.

Product Application Details	
Applications	Western Blot, Immunofluorescence
Recommended Dilutions	Western Blot 1:200, Immunofluorescence 1:50
Application Notes	WB: Detects a single 62 kDa protein representing nucleoporin p62 in rat liver nuclear envelope fractions. IF: Staining of p62 in rat liver with this antibody results in staining of both the nucleoplasmic and cytoplasmic surfaces of the nuclear pore complex (NPC).

Images

Western Blot: NUP62 Antibody (RL31) [NB120-2736] - Analysis of Nucleoporin p62 was performed by loading whole cell lysates (20ug) of the indicated samples and 5ul of prestained protein ladder onto a 4-20% Tris-Glycine polyacrylamide gel. Proteins were transferred to a nitrocellulose membrane, and blocked with 5% milk in TBST for 1 hour at room temperature. Nucleoporin p62 was detected at ~62kD in samples using a Nucleoporin p62 (RL31) monoclonal antibody at a dilution of 1:200 overnight at 4C on a rocking platform, followed by a goat anti-Armenian hamster IgG-HRP secondary antibody at a dilution of 1:10,000 for at least 30 minutes at room temperature.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB120-2736

NB120-5738	Goat anti-Armenian Hamster IgG (H+L) Secondary Antibody
NBP1-97046	Armenian Hamster IgG Isotype Control
H00023636-P01-10ug	Recombinant Human NUP62 GST (N-Term) Protein
NBP2-05223	NUP62 Overexpression Lysate

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB120-2736

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

