

# Product Datasheet

## Hepatitis C Virus Core Antigen Antibody (C7-50) [mFluor Violet 500 SE] NB120-2740MFV500

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NB120-2740MFV500](http://www.novusbio.com/NB120-2740MFV500)

Updated 9/20/2023 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NB120-2740MFV500](http://www.novusbio.com/reviews/destination/NB120-2740MFV500)



**NB120-2740MFV500**

Hepatitis C Virus Core Antigen Antibody (C7-50) [mFluor Violet 500 SE]

**Product Information**

<b>Unit Size</b>	0.1 ml
<b>Concentration</b>	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
<b>Storage</b>	Store at 4C in the dark.
<b>Clonality</b>	Monoclonal
<b>Clone</b>	C7-50
<b>Preservative</b>	0.05% Sodium Azide
<b>Isotype</b>	IgG1
<b>Conjugate</b>	mFluor Violet 500 SE
<b>Purity</b>	Immunogen affinity purified
<b>Buffer</b>	50mM Sodium Borate

**Product Description**


<b>Host</b>	Mouse
<b>Gene ID</b>	951475
<b>Gene Symbol</b>	POLY
<b>Species</b>	Human, Mouse, Virus
<b>Reactivity Notes</b>	Human reactivity reported in literature (PMID: 16940534). Mouse reactivity reported in literature (PMID: 22984515).
<b>Specificity/Sensitivity</b>	Detects hepatitis C virus (HCV) core protein from transfected human and primate cell lines.
<b>Immunogen</b>	Purified HCV core-GST fusion protein (genotype 1b).

**Product Application Details**

<b>Applications</b>	Western Blot, ELISA, Flow Cytometry, Immunohistochemistry, Immunoprecipitation, Immunofluorescence
<b>Recommended Dilutions</b>	Western Blot, Flow Cytometry, ELISA, Immunohistochemistry, Immunoprecipitation, Immunofluorescence
<b>Application Notes</b>	Optimal dilution of this antibody should be experimentally determined.

**Images**

Hepatitis C Virus Core Antigen Antibody (C7-50) [mFluor Violet 500 SE] - Vial of mFluor Violet 500 conjugated antibody. mFluor Violet 500 is optimally excited at 410 nm by the Violet laser (405 nm) and has an emission maximum of 501 nm.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NB120-2740MFV500](http://www.novusbio.com/reviews/submit/NB120-2740MFV500)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

