

# Product Datasheet

## PABP Antibody (10E10) [PerCP] NB120-6125PCP

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NB120-6125PCP](http://www.novusbio.com/NB120-6125PCP)

Updated 10/29/2025 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NB120-6125PCP](http://www.novusbio.com/reviews/destination/NB120-6125PCP)



**NB120-6125PCP**

PABP Antibody (10E10) [PerCP]

| Product Information         |                                                                                                                                                                                                                                                                            |
|-----------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size                   | 0.1 ml                                                                                                                                                                                                                                                                     |
| Concentration               | Please see the vial label for concentration. If unlisted please contact technical services.                                                                                                                                                                                |
| Storage                     | Store at 4C in the dark.                                                                                                                                                                                                                                                   |
| Clonality                   | Monoclonal                                                                                                                                                                                                                                                                 |
| Clone                       | 10E10                                                                                                                                                                                                                                                                      |
| Preservative                | 0.05% Sodium Azide                                                                                                                                                                                                                                                         |
| Isotype                     | IgG2b                                                                                                                                                                                                                                                                      |
| Conjugate                   | PerCP                                                                                                                                                                                                                                                                      |
| Purity                      | Protein G purified                                                                                                                                                                                                                                                         |
| Buffer                      | PBS                                                                                                                                                                                                                                                                        |
| Product Description         |                                                                                                                                                                                                                                                                            |
| Host                        | Mouse                                                                                                                                                                                                                                                                      |
| Gene ID                     | 26986                                                                                                                                                                                                                                                                      |
| Gene Symbol                 | PABPC1                                                                                                                                                                                                                                                                     |
| Species                     | Human, Bovine, Canine, Primate, Xenopus, Drosophila (Negative), Mouse (Negative)                                                                                                                                                                                           |
| Reactivity Notes            | Human, Monkey, Bovine, Dog, and Hamster and Xenopus laevis. Does not cross-react with Drosophila melanogaster or Mouse. Not yet tested in other species.                                                                                                                   |
| Immunogen                   | Recombinant human PABP protein. [Uniprot# P11940]                                                                                                                                                                                                                          |
| Product Application Details |                                                                                                                                                                                                                                                                            |
| Applications                | Western Blot, ELISA, Flow Cytometry, Immunocytochemistry/Immunofluorescence, Immunoprecipitation, Microarray, CyTOF-ready                                                                                                                                                  |
| Recommended Dilutions       | Western Blot, Flow Cytometry, ELISA, Immunocytochemistry/Immunofluorescence, Immunoprecipitation, Microarray, CyTOF-ready                                                                                                                                                  |
| Application Notes           | WB: Detects a band of approximately 70 kDa.<br><br>The observed molecular weight of the protein may vary from the listed predicted molecular weight due to post translational modifications, post translation cleavages, relative charges, and other experimental factors. |





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

[www.novusbio.com](http://www.novusbio.com)  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NB120-6125PCP**

---

|               |                                              |
|---------------|----------------------------------------------|
| NBP2-27231PCP | Mouse IgG2b Isotype Control (MPC-11) [PerCP] |
| NB120-6125UV  | PABP Antibody (10E10) [DyLight 350]          |
| NBP2-57869PEP | PABP Recombinant Protein Antigen             |
| NBP2-10839    | PABP Overexpression Lysate                   |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NB120-6125PCP](http://www.novusbio.com/reviews/submit/NB120-6125PCP)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

