

Product Datasheet

p53 Antibody [DyLight 350] NB200-171UV

Unit Size: 100 ul

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB200-171UV

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB200-171UV



NB200-171UV

p53 Antibody [DyLight 350]

| Product Information | |
|------------------------------------|--|
| Unit Size | 100 ul |
| Concentration | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage | Store at 4C in the dark. |
| Clonality | Polyclonal |
| Preservative | 0.05% Sodium Azide |
| Isotype | IgG |
| Conjugate | DyLight 350 |
| Purity | Immunogen affinity purified |
| Buffer | 50mM Sodium Borate |
| Product Description | |
| Description | This antibody can be used as the primary antibody in a PLA assay with the following as complementing antibodies: NB100-305, NB100-322, NBP1-19098, NBP1-19097, NB100-376, NB100-1803, NB100-224, NB100-61615, NB100-77315, NB500-160, NB100-381, NB100-381, NB100-596, NB200-347, NB100-97835, NBP1-06527, NBP1-065150 This conjugate is made on demand. Actual recovery may vary from the stated volume of this product. The volume will be greater than or equal to the unit size stated on the datasheet. |
| Host | Rabbit |
| Gene ID | 7157 |
| Gene Symbol | TP53 |
| Species | Human |
| Immunogen | The immunogen this antibody was made to, maps to a region between residues 50-100 of human tumor protein p53 using the SwissProt entry P04637 (GenelD 7157). |
| Notes | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries. |
| Product Application Details | |
| Applications | Immunohistochemistry, Chromatin Immunoprecipitation Sequencing |
| Recommended Dilutions | Immunohistochemistry, Chromatin Immunoprecipitation Sequencing |
| Application Notes | Optimal dilution of this antibody should be experimentally determined. |





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB200-171UV

| | |
|---------------|--|
| NBP2-24891UV | Rabbit IgG Isotype Control [DyLight 350] |
| NBP2-34495PEP | p53 Recombinant Protein Antigen |
| 1129-ER-050 | ErbB2/Her2 [Unconjugated] |
| DYC1043-2 | p53 [Biotin] |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB200-171UV

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

