

Product Datasheet

DDX5 Antibody NB200-351

Unit Size: 0.1 ml

Store at 4C. Do not freeze.

www.novusbio.com



technical@novusbio.com

Reviews: 1 Publications: 5

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB200-351

Updated 9/9/2025 v.20.1

**Earn rewards for product
reviews and publications.**

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB200-351



NB200-351**DDX5 Antibody**

Product Information	
Unit Size	0.1 ml
Concentration	0.2 mg/ml
Storage	Store at 4C. Do not freeze.
Clonality	Polyclonal
Preservative	0.09% Sodium Azide
Isotype	IgG
Purity	Immunogen affinity purified
Buffer	TBS and 0.1% BSA
Target Molecular Weight	69 kDa


Product Description	
Description	Novus Biologicals Rabbit DDX5 Antibody (NB200-351) is a polyclonal antibody validated for use in IHC, WB and IP. Anti-DDX5 Antibody: Cited in 5 publications. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	1655
Gene Symbol	DDX5
Species	Human, Mouse, Hamster
Reactivity Notes	Hamster reactivity reported in scientific literature (PMID: 28516729). Based on 100% sequence identity, this antibody is predicted to react with Orangutan and Chimpanzee
Immunogen	The immunogen recognized by this antibody maps to a region between residue 575 and the C-terminus (residue 614) of human DEAD (Asp-Glu-Ala-Asp) Polypeptide 5 using the numbering given in entry NP_004387.1 (GeneID 1655).

Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Immunohistochemistry, Immunoprecipitation
Recommended Dilutions	Western Blot 1:2500-1:10000, Immunohistochemistry 1:200 - 1:1000, Immunoprecipitation 2-5 ug/mg lysate, Immunohistochemistry-Paraffin 1 ug/ml
Application Notes	Optimal working dilutions should be determined experimentally by the investigator. Immunoprecipitation: the antibody may be diluted further than indicated. Epitope retrieval with citrate buffer pH6.0 is recommended for FFPE tissue sections.




Images


Western Blot: DDX5 Antibody [NB200-351] - Detection of Human and Mouse DDX5 by Western Blot. Samples: Whole cell lysate (50 ug) from HeLa, 293T, and mouse NIH3T3 cells prepared using NETN lysis buffer. Antibody: Affinity purified rabbit anti-DDX5 antibody NB200-351 used for WB at 0.1 ug/ml. Detection: Chemiluminescence with an exposure time of 1 second.




Immunohistochemistry: DDX5 Antibody [NB200-351] - Detection of human DDX5 by immunohistochemistry. Sample: FFPE section of human prostate carcinoma. Antibody: Affinity purified rabbit anti-DDX5 NB200-351 used at a dilution of 1:1000 (0.2ug/ml). Detection: DAB




Western Blot: DDX5 Antibody [NB200-351] - Detection of DDX5 in whole cell lysates from HeLa (5, 15, 50 ug), 293 (T; 50 ug), and MCF7 (M; 50 ug), and K562 (K; 50 ug) for WB (A) and HeLa (1 mg) for IP (B). Antibodies: NB 200-351 used at 0.04 ug/ml for WB (A & B) and at 3 ug/ml lysate for IP (B). DDX5 was also immunoprecipitated with NB 200-350 used at 3 ug/ml lysate. Detection: ECL 30 seconds.




Immunohistochemistry: DDX5 Antibody [NB200-351] - ~ ~Sample: FFPE section of human breast carcinoma.~ Antibody: Affinity purified rabbit anti-DDX5 ~Cat. No. Lot2 used at a dilution of 1:200 (1ug/ml). Detection: DAB. Counterstain: IHC Hematoxylin (blue).




Immunohistochemistry: DDX5 Antibody [NB200-351] - ~ Sample: FFPE section of mouse renal cell carcinoma.~ Antibody: Affinity purified rabbit anti-DDX5 ~Cat. No. Lot2 used at a dilution of 1:200 (1ug/ml). Detection: DAB. Counterstain: IHC Hematoxylin (blue).



Immunohistochemistry: DDX5 Antibody [NB200-351] - Detection of mouse DDX5 by immunohistochemistry. Sample: FFPE section of mouse skin. Antibody: Affinity purified rabbit anti-DDX5 NB200-351 used at a dilution of 1:1000 (0.2ug/ml). Detection: DAB



Immunoprecipitation: DDX5 Antibody [NB200-351] - Detection of human DDX5 by western blot of immunoprecipitates. Samples: Whole cell lysate (0.5 or 1.0 mg per IP reaction; 20% of IP loaded) from HeLa cells prepared using NETN lysis buffer. Antibodies: Affinity purified rabbit anti-DDX5 antibody NB200-351 (lot NB200-351-2) used for IP at 6 ug per reaction. DDX5 was also immunoprecipitated by a previous lot of this antibody (lot NB200-351-1) and another rabbit anti-DDX5 antibody. For blotting immunoprecipitated DDX5, NB200-351 was used at 1 ug/ml. Detection: Chemiluminescence with an exposure time of 1 second.



Publications

Sun X, Wang S, Lin X et al. Proteome analysis of Duck Tembusu virus (DTMUV)-infected BHK-21 cells. *Proteomics* 2017-06-20 [PMID: 28516729] (WB, Hamster)

Muller RY, Hammond MC, Rio DC, Lee YJ. An Efficient Method for Electroporation of Small Interfering RNAs into ENCODE Project Tier 1 GM12878 and K562 Cell Lines. *J Biomol Tech.* 2015-10-29 [PMID: 26543439] (WB, Human)

Saporita AJ, Chang HC, Winkeler CL et al. RNA helicase DDX5 is a p53-independent target of ARF that participates in ribosome biogenesis. *Cancer Res* 2011-11-01 [PMID: 21937682]

Guo J, Hong F, Loke J et al. A DDX5 S480A polymorphism is associated with increased transcription of fibrogenic genes in hepatic stellate cells. *J Biol Chem* 2010-02-01 [PMID: 20022962]

Chiu YL, Witkowska HE, Hall SC et al. High-molecular-mass APOBEC3G complexes restrict Alu retrotransposition. *Proc Natl Acad Sci U S A* 2006-10-01 [PMID: 17030807]



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB200-351

NB800-PC1	HeLa Whole Cell Lysate
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB200-351

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

