

Product Datasheet

Peripherin Antibody (8G2) [DyLight 350] NB300-138UV-0.1ml

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB300-138UV

Updated 10/30/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB300-138UV



NB300-138UV-0.1ml

Peripherin Antibody (8G2) [DyLight 350]

Product Information	
Unit Size	0.1 ml
Concentration	This product is unpurified. The exact concentration of antibody is not quantifiable.
Storage	Store at 4C in the dark.
Clonality	Monoclonal
Clone	8G2
Preservative	0.05% Sodium Azide
Isotype	IgG1
Conjugate	DyLight 350
Purity	Tissue culture supernatant
Buffer	50mM Sodium Borate
Target Molecular Weight	57 kDa
Product Description	
Host	Mouse
Gene ID	5630
Gene Symbol	PRPH
Species	Human, Mouse, Rat, Porcine, Bovine
Marker	Un-myelinated Neurons Marker
Immunogen	Recombinant rat Peripherin purified from E. coli. [UniProt# P21807]
Notes	Dylight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries. This conjugate is made on demand. Actual recovery may vary from the stated volume of this product. The volume will be greater than or equal to the unit size stated on the datasheet.
Product Application Details	
Applications	Western Blot, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry
Recommended Dilutions	Western Blot, Immunohistochemistry, Immunocytochemistry/ Immunofluorescence
Application Notes	This Peripherin antibody can be used for Immunocytochemistry/Immunofluorescence, ELISA and Western blotting. In WB a band can be seen at ~57 kDa. The observed molecular weight of the protein may vary from the listed predicted molecular weight due to post translational modifications, post translation cleavages, relative charges, and other experimental factors.

**Novus Biologicals USA**

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB300-138UV-0.1ml

NBP1-97005UV	Mouse IgG1 Isotype Control (MG1) [DyLight 350]
NB300-138FR-0.1ml	Peripherin Antibody (8G2) [DyLight 680]
H00005630-P01-10ug	Recombinant Human Peripherin GST (N-Term) Protein
256-GF-100	beta-NGF [Unconjugated]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB300-138UV

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

