

Product Datasheet

Endothelin-1 Antibody (TR.ET.48.5) [Allophycocyanin] NB300-526APC

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB300-526APC

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB300-526APC



NB300-526APC

Endothelin-1 Antibody (TR.ET.48.5) [Allophycocyanin]

| Product Information | |
|------------------------------------|--|
| Unit Size | 0.1 ml |
| Concentration | Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services. |
| Storage | Store at 4C in the dark. |
| Clonality | Monoclonal |
| Clone | TR.ET.48.5 |
| Preservative | 0.05% Sodium Azide |
| Isotype | IgG1 |
| Conjugate | Allophycocyanin |
| Purity | Protein A or G purified |
| Buffer | PBS |
| Product Description | |
| Host | Mouse |
| Gene ID | 1906 |
| Gene Symbol | EDN1 |
| Species | Human, Mouse, Rat, Porcine, Canine, Rabbit, Sheep |
| Reactivity Notes | Canine reactivity reported in scientific literature (PMID: 11673223). Rabbit reactivity reported in scientific literature (PMID: 12181591). Porcine reactivity reported in scientific literature (PMID: 21587116). Please note that this antibody is reactive to Mouse and derived from the same host, Mouse. Additional Mouse on Mouse blocking steps may be required for IHC and ICC experiments. Please contact Technical Support for more information. |
| Immunogen | Pure endothelin-1 conjugated to KLH. |
| Product Application Details | |
| Applications | Western Blot, Flow Cytometry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Frozen, Immunohistochemistry-Paraffin, Radioimmunoassay |
| Recommended Dilutions | Western Blot, Flow Cytometry, Immunohistochemistry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry-Paraffin, Immunohistochemistry-Frozen, Radioimmunoassay |





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB300-526APC

| | |
|--------------------|---|
| NBP1-97005APC | Mouse IgG1 Isotype Control (MG1) [Allophycocyanin] |
| NB300-526UV | Endothelin-1 Antibody (TR.ET.48.5) [DyLight 350] |
| H00001906-Q01-10ug | Recombinant Human Endothelin-1 GST (N-Term) Protein |
| 210-TA-005 | TNF-alpha [Unconjugated] |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB300-526APC

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

