

Product Datasheet

Dkk-2 Antibody NB300-810

Unit Size: 0.1 mg

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB300-810

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB300-810



NB300-810**Dkk-2 Antibody****Product Information**

| | |
|----------------------|--|
| Unit Size | 0.1 mg |
| Concentration | 0.5 mg/ml |
| Storage | Store at -20C. Avoid freeze-thaw cycles. |
| Clonality | Polyclonal |
| Preservative | 0.02% Sodium Azide |
| Isotype | IgG |
| Purity | Immunogen affinity purified |
| Buffer | Tris saline (20 mM Tris pH 7.3, 150 mM NaCl), 0.5% BSA |

Product Description


| | |
|--------------------------------|--|
| Description | Novus Biologicals Goat Dkk-2 Antibody (NB300-810) is a polyclonal antibody validated for use in WB and ELISA. Anti-Dkk-2 Antibody: Cited in 1 publication. All Novus Biologicals antibodies are covered by our 100% guarantee. |
| Host | Goat |
| Gene ID | 27123 |
| Gene Symbol | DKK2 |
| Species | Human |
| Specificity/Sensitivity | From sequence analysis it is possible that this product will crossreact with DKK4, unfortunately this has not yet been specifically tested. |
| Immunogen | Peptide with sequence ATYSSKARLHVCQKI corresponding to C-Terminus according to NP_055236.1. |

Product Application Details

| | |
|------------------------------|---|
| Applications | Western Blot, Peptide ELISA |
| Recommended Dilutions | Western Blot 1 - 3 ug/ml, Peptide ELISA Detection limit 1:8000 |
| Application Notes | WB: Approx. 28 kDa band observed in human heart lysate (calculated MW of 28.5 kDa band according to NP_055236). |

Images

Western Blot: Dkk-2 Antibody [NB300-810] - NB300-810 (2 ug/ml) staining of human heart lysate (35 ug protein in RIPA buffer). Primary incubation was 1 hour. Detected by chemiluminescence.

**Publications**

Krupnik VE, Sharp JD, Jiang C et al. Functional and structural diversity of the human Dickkopf gene family. Gene 1999-10-01 [PMID: 10570958]



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB300-810

| | |
|-----------------|--|
| HAF017 | Rabbit anti-Goat IgG Secondary Antibody [HRP (Horseradish Peroxidase)] |
| HAF109 | Donkey anti-Goat IgG Secondary Antibody [HRP (Horseradish Peroxidase)] |
| NB410-28088-1mg | Goat IgG Isotype Control |
| NBP2-68703PEP | Dkk-2 Recombinant Protein Antigen |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB300-810

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications



