

Product Datasheet

TXNDC Antibody NB300-995

Unit Size: 0.1 mg

Store at -20C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB300-995

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB300-995



NB300-995

TXNDC Antibody

| Product Information | |
|---------------------|--|
| Unit Size | 0.1 mg |
| Concentration | 0.5 mg/ml |
| Storage | Store at -20C. Avoid freeze-thaw cycles. |
| Clonality | Polyclonal |
| Preservative | 0.02% Sodium Azide |
| Isotype | IgG |
| Purity | Immunogen affinity purified |
| Buffer | Tris saline (20 mM Tris pH 7.3, 150 mM NaCl), 0.5% BSA |

| Product Description | |
|---------------------|--|
| Host | Goat |
| Gene ID | 81542 |
| Gene Symbol | TMX1 |
| Species | Human |
| Immunogen | Peptide with sequence C-RSLGPSLATDKS corresponding to C-Terminus according to NP_110382.3. |

| Product Application Details | |
|-----------------------------|--|
| Applications | Western Blot, Peptide ELISA |
| Recommended Dilutions | Western Blot 0.01 - 0.03 ug/mL, Peptide ELISA Detection limit 1:64000 |
| Application Notes | WB: Approx. 35 kDa band observed in human liver lysates (calculated MW of 31.8 kDa band according to NP_110382). IHC-P: Human spleen shows textured staining in the cytoplasm of some splenocytes. |

Images

Western Blot: TXNDC Antibody [NB300-995] - Staining of Human Liver lysate (35 ug protein in RIPA buffer). Antibody at 0.01 ug/mL. Detected by chemiluminescence.

**Publications**

Matsuo Y, Akiyama N, Nakamura H et al. Identification of a novel thioredoxin-related transmembrane protein. J Biol Chem 2001-03-30 [PMID: 11152479]



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB300-995

| | |
|-----------------|--|
| NBL1-17463 | TXNDC Overexpression Lysate |
| HAF017 | Rabbit anti-Goat IgG Secondary Antibody [HRP (Horseradish Peroxidase)] |
| HAF109 | Donkey anti-Goat IgG Secondary Antibody [HRP (Horseradish Peroxidase)] |
| NB410-28088-1mg | Goat IgG Isotype Control |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB300-995

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications



