

Product Datasheet

Goat IgG Isotype Control [DyLight 650] NB410-28088C-500ul

Unit Size: 500 ul

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB410-28088C

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB410-28088C



NB410-28088C-500ul

Goat IgG Isotype Control [DyLight 650]

Product Information	
Unit Size	500 ul
Concentration	Please contact technical services for concentration.
Storage	Store at 4C in the dark.
Clonality	Polyclonal
Preservative	0.05% Sodium Azide
Isotype	IgG
Conjugate	DyLight 650
Buffer	50 mM Sodium Borate.

Product Description	
Description	This goat IgG immunoglobulin is useful as an isotype control for goat IgG antibodies.
Host	Goat
Gene ID	3500
Gene Symbol	IGHG1
Species	Goat
Notes	DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.

Product Application Details	
Applications	ELISA, Fluorophore-linked immunosorbent assay, Flow Cytometry
Recommended Dilutions	Flow Cytometry 1:10 - 1:1000, ELISA 1:100 - 1:2000, Fluorophore-linked immunosorbent assay 1:100 - 1:2000
Application Notes	Optimal dilution of this antibody should be experimentally determined. Use the isotype control at the same concentration as the experimental antibody. Optimal dilution of this antibody should be experimentally determined.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB410-28088C-500ul

NBP1-99014-100ug	Recombinant Mouse IgG His Protein
NB200-540	Complement C3 Antibody (11H9) - BSA Free
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NBP2-29373-100Tests	Annexin V Apoptosis Kit [FITC]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Isotype Controls are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB410-28088C

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

