

# Product Datasheet

## Goat IgG Isotype Control [PerCP] NB410-28088PCP

Unit Size: 1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NB410-28088PCP](http://www.novusbio.com/NB410-28088PCP)

Updated 10/23/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NB410-28088PCP](http://www.novusbio.com/reviews/destination/NB410-28088PCP)



**NB410-28088PCP**

Goat IgG Isotype Control [PerCP]

<b>Product Information</b>	
<b>Unit Size</b>	1 ml
<b>Concentration</b>	Please contact technical services for concentration.
<b>Storage</b>	Store at 4C in the dark.
<b>Clonality</b>	Polyclonal
<b>Preservative</b>	0.05% Sodium Azide
<b>Isotype</b>	IgG
<b>Conjugate</b>	PerCP
<b>Buffer</b>	PBS (Low Endotoxin).

<b>Product Description</b>	
<b>Description</b>	This goat IgG immunoglobulin is useful as an isotype control for goat IgG antibodies.
<b>Host</b>	Goat
<b>Gene ID</b>	3500
<b>Gene Symbol</b>	IGHG1
<b>Species</b>	Goat

<b>Product Application Details</b>	
<b>Applications</b>	ELISA, Fluorophore-linked immunosorbent assay, Flow Cytometry
<b>Recommended Dilutions</b>	Flow Cytometry 1:10 - 1:1000, ELISA 1:100 - 1:2000, Fluorophore-linked immunosorbent assay 1:100 - 1:2000
<b>Application Notes</b>	Optimal dilution of this antibody should be experimentally determined. Use the isotype control at the same concentration as the experimental antibody. Optimal dilution of this antibody should be experimentally determined.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NB410-28088PCP**

---

NBP1-99013-50ug	Recombinant Mouse IgG Protein
NB200-540	Complement C3 Antibody (11H9) - BSA Free
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NBP2-29373-100Tests	Annexin V Apoptosis Kit [FITC]

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Isotype Controls are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NB410-28088PCP](http://www.novusbio.com/reviews/submit/NB410-28088PCP)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

