

Product Datasheet

UBE2R1/CDC34 Antibody (14B5) - Azide and BSA Free NB500-168

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB500-168

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB500-168



NB500-168

UBE2R1/CDC34 Antibody (14B5) - Azide and BSA Free

Product Information	
Unit Size	0.1 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	14B5
Preservative	No Preservative
Isotype	IgG1
Purity	Protein G purified
Buffer	PBS (pH 7.2)
Target Molecular Weight	27 kDa


Product Description	
Description	Novus Biologicals Mouse UBE2R1/CDC34 Antibody (14B5) - Azide and BSA Free (NB500-168) is a monoclonal antibody validated for use in IHC, WB, ICC/IF and IP. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	997
Gene Symbol	CDC34
Species	Human
Immunogen	Full-length human UBE2R1/CDC34 expressed in E. coli.

Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Immunocytochemistry/Immunofluorescence, Immunohistochemistry, Immunoprecipitation
Recommended Dilutions	Western Blot, Immunohistochemistry 10 - 1:500, Immunocytochemistry/Immunofluorescence, Immunoprecipitation, Immunohistochemistry-Paraffin



Images


Western Blot: UbchH3/Cdc34 Antibody (14B5) [NB500-168] - Sample (30 ug of whole cell lysate) A: A549 B: HCT116 10% SDS PAGE; antibody diluted at 1:1000.



Immunocytochemistry/Immunofluorescence: UbchH3/Cdc34 Antibody (14B5) [NB500-168] - Analysis of paraformaldehyde-fixed HeLa, using CDC34 antibody at 1:100 dilution.



Immunohistochemistry-Paraffin: UbchH3/Cdc34 Antibody (14B5) [NB500-168] - Analysis of paraffin-embedded HBL435 xenograft, using CDC34 antibody at 1:200 dilution.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB500-168

NB800-PC6	293 Whole Cell Lysate
HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-97005-0.5mg	Mouse IgG1 Isotype Control (MG1)

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB500-168

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

