

Product Datasheet

Survivin Antibody (32.1) [DyLight 550] NB500-237R

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB500-237R

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB500-237R



NB500-237R

Survivin Antibody (32.1) [DyLight 550]

| Product Information | |
|------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Size | 0.1 ml |
| Concentration | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage | Store at 4C in the dark. |
| Clonality | Monoclonal |
| Clone | 32.1 |
| Preservative | 0.05% Sodium Azide |
| Isotype | IgG1 Kappa |
| Conjugate | DyLight 550 |
| Purity | Protein G purified |
| Buffer | 50mM Sodium Borate |
| Target Molecular Weight | 16 kDa |
| Product Description | |
| Host | Mouse |
| Gene ID | 332 |
| Gene Symbol | BIRC5 |
| Species | Human, Hamster |
| Reactivity Notes | Hamster reactivity reported in scientific literature (PMID: 23405201). |
| Specificity/Sensitivity | Survivin Antibody (32.1) [NB500-237] is specific for the nuclear form of survivin. |
| Immunogen | This Survivin Antibody (32.1) [DyLight 550] was developed against full length recombinant human Survivin [UniProt# O15392] |
| Notes | DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries. |
| Product Application Details | |
| Applications | Western Blot, ELISA, Flow Cytometry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin, Immunoprecipitation |
| Recommended Dilutions | Western Blot, Flow Cytometry, ELISA, Immunohistochemistry, Immunocytochemistry/ Immunofluorescence, Immunoprecipitation, Immunohistochemistry-Paraffin |
| Application Notes | Optimal dilution of this antibody should be experimentally determined. |





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB500-237R

| | |
|----------------|-------------------------------------------------------------|
| NBP1-43319R | Mouse IgG1 Kappa Isotype Control (P3.6.2.8.1) [DyLight 550] |
| NB500-237H | Survivin Antibody (32.1) [HRP] |
| NB110-92717PEP | Survivin Antibody Blocking Peptide |
| DVE00 | VEGF [HRP] |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB500-237R

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

