

Product Datasheet

Survivin Antibody (60.11) [DyLight 650] NB500-238C

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB500-238C

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB500-238C



NB500-238C

Survivin Antibody (60.11) [DyLight 650]

Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark.
Clonality	Monoclonal
Clone	60.11
Preservative	0.05% Sodium Azide
Isotype	IgG2a Kappa
Conjugate	DyLight 650
Purity	Protein G purified
Buffer	50mM Sodium Borate
Target Molecular Weight	16 kDa

Product Description	
Host	Mouse
Gene ID	332
Gene Symbol	BIRC5
Species	Human, Mouse, Rat
Specificity/Sensitivity	Survivin Antibody (60.11) [NB500-238] is specific for the cytoplasmic form of Survivin.
Immunogen	This Survivin Antibody (60.11) [DyLight 650] was developed against full length recombinant human Survivin [UniProt# O15392]
Notes	DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.

Product Application Details	
Applications	Western Blot, ELISA, Flow Cytometry, Flow (Intracellular), Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Frozen, Immunohistochemistry-Paraffin, Proximity Ligation Assay, Knockdown Validated
Recommended Dilutions	Western Blot, Flow Cytometry, ELISA, Immunohistochemistry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry-Paraffin, Immunohistochemistry-Frozen, Proximity Ligation Assay, Flow (Intracellular), Knockdown Validated
Application Notes	Optimal dilution of this antibody should be experimentally determined.

Publications

Kusner LL, Ciesielski MJ, Marx A et al. Survivin as a potential mediator to support autoreactive cell survival in myasthenia gravis: a human and animal model study. PLoS ONE 2014-07-23 [PMID: 25050620] (FLOW, Human)

Details:

Using the DyLight 650 conjugated version of NB500-238, catalog number NB500-238C.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB500-238C

NBP1-96981C	Mouse IgG2a Kappa Isotype Control (M2AK) [DyLight 650]
NB500-238F	Survivin Antibody (60.11) [FITC]
NB110-92717PEP	Survivin Antibody Blocking Peptide
DVE00	VEGF [HRP]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB500-238C

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

