

# Product Datasheet

## L-Selectin/CD62L Antibody (IVA94) [FITC] NB500-497F

Unit Size: 0.1 ml

Store at 4C in the dark.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NB500-497F](http://www.novusbio.com/NB500-497F)

Updated 10/23/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NB500-497F](http://www.novusbio.com/reviews/destination/NB500-497F)



**NB500-497F****L-Selectin/CD62L Antibody (IVA94) [FITC]**

Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark.
Clonality	Monoclonal
Clone	IVA94
Preservative	0.05% Sodium Azide
Isotype	IgG1
Conjugate	FITC
Purity	Protein A or G purified
Buffer	PBS

Product Description	
Host	Mouse
Gene ID	6402
Gene Symbol	SELL
Species	Human, Bovine
Specificity/Sensitivity	The antibody IVA94 reacts with CD62L antigen (bovine). CD62L (L-selectin) is a 74-95 kDa single chain type I glycoprotein expressed on most peripheral blood B lymphocytes, T lymphocytes, monocytes and granulocytes; it is also present on a subset of NK cells and certain hematopoietic malignant cells.
Immunogen	Bovine leukocytes.

Product Application Details	
Applications	Flow Cytometry, Immunoprecipitation, CyTOF-ready, Immunocytochemistry
Recommended Dilutions	Flow Cytometry, Immunoprecipitation, Immunocytochemistry, CyTOF-ready
Application Notes	Optimal dilution of this antibody should be experimentally determined.

**Images**

L-Selectin/CD62L Antibody (IVA94) [FITC] [NB500-497F] - Vial of FITC conjugated antibody. FITC is optimally excited at 498 nm by the Blue laser (488 nm) and has an emission maximum of 519 nm.



FITC

LASER (nm)	FILTER
Blue (488)	525/50
EXCITATION MAX (nm)	EMISSION MAX (nm)
498	519



### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NB500-497F**

---

NBP1-96849	Mouse IgG1 Isotype Control (MG1) [FITC]
NBP2-76545PEP	L-Selectin/CD62L Recombinant Protein Antigen
210-TA-005	TNF-alpha [Unconjugated]
728-LS-100	L-Selectin/CD62L

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NB500-497F](http://www.novusbio.com/reviews/submit/NB500-497F)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

