

Product Datasheet

HDAC10 Antibody - BSA Free

NB600-223

Unit Size: 0.1 mg

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Reviews: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB600-223

Updated 4/13/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB600-223



NB600-223


HDAC10 Antibody - BSA Free

Product Information	
Unit Size	0.1 mg
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.02% Sodium Azide
Isotype	IgG
Purity	Immunogen affinity purified
Buffer	PBS
Target Molecular Weight	62 kDa
Product Description	
Host	Rabbit
Gene ID	83933
Gene Symbol	HDAC10
Species	Human
Reactivity Notes	Human.
Specificity/Sensitivity	Reacts with the human HDAC 10 protein.
Immunogen	Synthetic peptide corresponding to a portion of human HDAC 10.
Product Application Details	
Applications	Western Blot, Simple Western
Recommended Dilutions	Western Blot 1:500-1:1000, Simple Western
Application Notes	<p>Western blot, a predominant band is seen at ~65-70 kDa; it has been reported that HDAC10 variants seem present in certain cell lines and tissues that are as of yet uncharacterized. Other bands may be seen below 50 kDa, and in some cases they are predominant. Simple Western reported in verified customer review.</p> <p>See Simple Western Antibody Database for Simple Western validation: Tested in Human SVG-A Protein Lysate, separated by Size, antibody dilution of 1:25</p>




Images

Western Blot: HDAC10 Antibody [NB600-223] - WB analysis of HDAC10 in HeLa cells with NB600-223 at a dilution of 1:250.



Simple Western: HDAC10 Antibody [NB600-223] - 1:25 dilution against the indicated concentration of protein lysate (Human SVG-A). Image from verified customer review.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB600-223

NBL1-11478	HDAC10 Overexpression Lysate
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB600-223

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

