

Product Datasheet

Aggrecan Antibody (BC-3) [DyLight 550] NB600-504R-100ul

Unit Size: 100 ul

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB600-504R

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB600-504R



NB600-504R-100ul

Aggrecan Antibody (BC-3) [DyLight 550]

Product Information	
Unit Size	100 ul
Concentration	This product is unpurified. The exact concentration of antibody is not quantifiable.
Storage	Store at 4C in the dark.
Clonality	Monoclonal
Clone	BC-3
Preservative	0.05% Sodium Azide
Isotype	IgG1 Kappa
Conjugate	DyLight 550
Purity	Tissue culture supernatant
Buffer	50mM Sodium Borate
Target Molecular Weight	250 kDa
Product Description	
Host	Mouse
Gene ID	176
Gene Symbol	ACAN
Species	Human, Mouse, Rat, Porcine, Bovine, Canine, Equine, Feline, Guinea Pig, Rabbit, Sheep
Reactivity Notes	Cross-reacts with Human, Cow, Cat, Dog, Guinea pig, Horse, Pig, Rabbit, Rat and Sheep. Mouse reactivity reported in scientific literature (PMID: 27185069). Not yet tested in other species.
Marker	Chondrocyte Marker
Specificity/Sensitivity	Recognizes the aggrecanase (ADAMTS-1, -4 & -5)-generated N-terminal neopeptide ARG after cleavage between amino acids EGE and ARG within the interglobular domain of aggrecanase-catabolised aggrecan.
Immunogen	ARGSV... synthetic peptide conjugate.
Notes	DyLight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.
Product Application Details	
Applications	Western Blot, ELISA, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Frozen, Immunohistochemistry-Paraffin
Recommended Dilutions	Western Blot, ELISA, Immunohistochemistry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry-Paraffin, Immunohistochemistry-Frozen
Application Notes	This antibody should work in IHC on formalin- or paraformaldehyde-fixed paraffin embedded sections as well as either alcohol-fixed frozen sections or un-fixed snap-frozen sections. This antibody can also be used with ELISA and Western blot, where it detects a variety of epitopes between 50 and 250 kDa. The observed molecular weight of the protein may vary from the listed predicted molecular weight due to post translational modifications, post translation cleavages, relative charges, and other experimental factors.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB600-504R-100ul

NBP1-43319R	Mouse IgG1 Kappa Isotype Control (P3.6.2.8.1) [DyLight 550]
NB600-504B-100ul	Aggrecan Antibody (BC-3) [Biotin]
NBP2-55540PEP	Aggrecan Recombinant Protein Antigen
210-TA-005	TNF-alpha [Unconjugated]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB600-504R

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

