

Product Datasheet

alpha Tubulin [ac Lys40] Antibody (6-11B-1) NB600-567

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Reviews: 1 Publications: 2


Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB600-567

Updated 11/11/2025 v.20.1

**Earn rewards for product
reviews and publications.**

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB600-567



NB600-567

alpha Tubulin [ac Lys40] Antibody (6-11B-1)

Product Information	
Unit Size	0.1 ml
Concentration	This product is unpurified. The exact concentration of antibody is not quantifiable.
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	6-11B-1
Preservative	0.09% Sodium Azide
Isotype	IgG2b
Purity	Unpurified
Buffer	Ascites
Target Molecular Weight	50 kDa
Product Description	
Host	Mouse
Gene ID	7846
Gene Symbol	TUBA1A
Species	Human, Mouse, Rat, Porcine, Amphibian, Bovine, Canine, Chicken, Hamster, Invertebrate, Parasite, Plant, Primate, Protozoa
Reactivity Notes	Please note that this antibody is reactive to Mouse and derived from the same host, Mouse. Additional Mouse on Mouse blocking steps may be required for IHC and ICC experiments. Please contact Technical Support for more information.
Marker	Microtubule Marker
Specificity/Sensitivity	This alpha Tubulin [ac Lys40] Antibody (6-11B-1) recognizes an epitope located on the alpha 3 isoform of Chlamydomonas axonemal alpha-tubulin, within four residues of Lys-40 when this amino acid is acetylated
Immunogen	This alpha Tubulin [ac Lys40] Antibody (6-11B-1) was developed against acetylated tubulin from the outer arm of Strongylocentrotus purpuratus (sea urchin)
Product Application Details	
Applications	Western Blot, Dot Blot, ELISA, Electron Microscopy, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunoprecipitation, Radioimmunoassay
Recommended Dilutions	Western Blot 1:2000, ELISA 1:100-1:2000, Immunohistochemistry 1:10-1:500, Immunocytochemistry/ Immunofluorescence 1:10-1:2000, Immunoprecipitation, Radioimmunoassay 1:100-1:2000, Dot Blot 1:100-1:2000, Electron Microscopy 1:10-1:500




Images


Immunohistochemistry: alpha Tubulin [ac Lys40] Antibody (6-11B-1) [NB600-567] - Nerve pathways in the embryonic zebrafish brain and spinal cord revealed by alpha Tubulin (ac Lys40) Antibody (6-11B-1). Immunoperoxidase staining of a whole-mount embryo at 24 hours of development.




Western Blot: alpha Tubulin [ac Lys40] Antibody (6-11B-1) [NB600-567] - Western Blot validation of monoclonal anti-Acetylated Tubulin antibody clone: 6-11B-1 at 1:2000 dilution with goat anti-Mouse IgG-Peroxidase secondary antibody/chemiluminescent substrate development. The lysates employed were from cell lines: 1. HeLa, 2. U87, 3. COS7, 4. P19, 5. RAT2, 6. CHO, 7. MDBK and 8. MDCK Alpha tubulin molecular weight: 50 kDa.



Immunocytochemistry/Immunofluorescence: alpha Tubulin [ac Lys40] Antibody (6-11B-1) [NB600-567] - HeLa cells were fixed with 4% paraformaldehyde followed permeabilization with 0.5% Triton X-100. Fixed cells were stained with 1:2,000 monoclonal Tubulin, Acetylated antibody Clone: 6-11B-1. The antibody was detected using goat anti-Mouse IgG, Cy3 conjugate followed by counterstaining with DAPI (blue).



Immunoprecipitation: alpha Tubulin [ac Lys40] Antibody (6-11B-1) [NB600-567] - Identification of a surrogate peptide YCLEHGIQPDGQMPSDK40(ac)TIGGGDDSF with mass spectrometry. This peptide was immunoprecipitated with 6-11B-1 mAb from peptic digested rat brain tissue lysate. Image submitted via customer review.



Publications

Mardilovich K, Gabrielsen M et al. Elevated LIM kinase 1 in nonmetastatic prostate cancer reflects its role in facilitating androgen receptor nuclear translocation. Mol Cancer Ther 2015-01-01 [PMID: 25344584] (ICC/IF, Human)

Yang X, Naughton S, Ha Z et al. Mass Spectrometric Quantitation of Tubulin Acetylation from Pepsin-Digested Rat Brain Tissue Using a Novel Stable Isotope Standard and Capture by Anti-Peptide Antibody (SISCAPA) Method. Anal Chem 2018-01-24 [PMID: 29320166] (IP, Rat)



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB600-567

Earn gift cards/discounts by submitting a publication using this product:

www.novusbio.com/publications

