

# Product Datasheet

## Goat anti-Feline IgM Heavy Chain Secondary Antibody [FITC] NB7259

Unit Size: 1 ml

Store at 4C. Do not freeze.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

**Reviews: 1**


Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NB7259](http://www.novusbio.com/NB7259)

Updated 9/22/2021 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NB7259](http://www.novusbio.com/reviews/destination/NB7259)



**NB7259****Goat anti-Feline IgM Heavy Chain Secondary Antibody [FITC]**


| Product Information |                                                    |
|---------------------|----------------------------------------------------|
| Unit Size           | 1 ml                                               |
| Concentration       | 1.0 mg/ml                                          |
| Storage             | Store at 4C. Do not freeze.                        |
| Clonality           | Polyclonal                                         |
| Preservative        | 0.09% Sodium Azide                                 |
| Isotype             | IgG                                                |
| Conjugate           | FITC                                               |
| Purity              | Immunogen affinity purified                        |
| Buffer              | Phosphate Buffered Saline (PBS) containing 0.2%BSA |

| Product Description     |                                                                                                                             |
|-------------------------|-----------------------------------------------------------------------------------------------------------------------------|
| Host                    | Goat                                                                                                                        |
| Gene ID                 | 3507                                                                                                                        |
| Species                 | Feline                                                                                                                      |
| Specificity/Sensitivity | By immunoelectrophoresis and ELISA this reacts specifically with cat IgM. This may cross react with IgM from other species. |
| Immunogen               | Cat IgM                                                                                                                     |

| Product Application Details |                                                                                                                                                    |
|-----------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------|
| Applications                | Flow Cytometry, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin                                       |
| Recommended Dilutions       | Flow Cytometry 1:10-1:1000, Immunohistochemistry 1:50 - 1:500, Immunocytochemistry/ Immunofluorescence 1:50 - 1:500, Immunohistochemistry-Paraffin |

**Images**

Flow Cytometry: Goat anti-Feline IgM Heavy Chain Secondary Antibody [FITC] [NB7259] - The cat cells were washed and Incubated with the IgM antibody. It was used to detect Bcells. This image was submitted via customer review.





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NB7259**

---

|             |                                       |
|-------------|---------------------------------------|
| NB100-524   | NOD2 Antibody (2D9) - BSA Free        |
| NBP1-96772  | Human IgM Heavy Chain Isotype Control |
| 7268-CT-100 | CTLA-4 [Unconjugated]                 |
| NB200-541   | Complement C4b/d Antibody (16D2)      |

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Secondary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NB7259](http://www.novusbio.com/reviews/submit/NB7259)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

