

Product Datasheet

Goat anti-Human Kappa Light Chain Secondary Antibody NB7480

Unit Size: 1 ml

Store at 4C. Do not freeze.

www.novusbio.com



technical@novusbio.com

Publications: 2


Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB7480

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB7480



NB7480**Goat anti-Human Kappa Light Chain Secondary Antibody**

Product Information	
Unit Size	1 ml
Concentration	1.0 mg/ml
Storage	Store at 4C. Do not freeze.
Clonality	Polyclonal
Preservative	0.09% Sodium Azide
Isotype	IgG
Purity	Immunogen affinity purified
Buffer	Phosphate Buffered Saline (PBS)

Product Description	
Host	Goat
Gene ID	3514
Gene Symbol	IGKC
Species	Human
Specificity/Sensitivity	By immunoelectrophoresis and ELISA this reacts specifically with kappa free light chains. No was detected against lambda light chain or other serum proteins.
Immunogen	Human Kappa Free Light Chain

Product Application Details	
Applications	Western Blot, ELISA, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin
Recommended Dilutions	Western Blot 1:1000-1:30000, ELISA 1:1000-1:30000 (Primary) 1:100-1:500 (for coating plates), Immunohistochemistry 1:200- 1:2000, Immunocytochemistry/ Immunofluorescence 1:200- 1:2000, Immunohistochemistry-Paraffin

Publications

Rocchetti MT, Papale M, d'Apollo AM et al. Association of urinary laminin G-like 3 and free K light chains with disease activity and histological injury in IgA nephropathy Clin J Am Soc Nephrol 2013-04-20 [PMID: 23599406]

Ishii Y, Tsukahara M, Wakamatsu K. A rapid method for simultaneous evaluation of free light chain content and aggregate content in culture media of Chinese hamster ovary cells expressing monoclonal antibodies for cell line screening. J Biosci Bioeng 2016-04-01 [PMID: 26467692]





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NB7480

H00003514-P01-10ug	Recombinant Human Kappa Light Chain GST (N-Term) Protein
NB100-1028	AIF-1/Iba1 Antibody
NBP2-34258-0.1mg	Kappa Light Chain Antibody (HP6053 + L1C1)
243-B3-002	TGF-beta 3 [Unconjugated]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Secondary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB7480

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

