

Product Datasheet

PD-1 Antibody (4F12) [Allophycocyanin/Cy7] NBP1-04328APCCY7

Unit Size: 0.1 ml

Store at 4C in the dark. Do not freeze.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-04328APCCY7

Updated 10/23/2024 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-04328APCCY7



NBP1-04328APCCY7

PD-1 Antibody (4F12) [Allophycocyanin/Cy7]

Product Information	
Unit Size	0.1 ml
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark. Do not freeze.
Clonality	Monoclonal
Clone	4F12
Preservative	0.05% Sodium Azide
Isotype	IgG2b Kappa
Conjugate	Allophycocyanin/Cy7
Purity	Protein G purified
Buffer	PBS

Product Description	
Host	Mouse
Gene ID	5133
Gene Symbol	PDCD1
Species	Human, Canine
Reactivity Notes	Use in Canine reported in scientific literature (PMID:34923192) .
Immunogen	Recombinant human PD-1 (21-167aa) purified from E. coli

Product Application Details	
Applications	Western Blot, ELISA, Flow Cytometry, Immunocytochemistry/ Immunofluorescence
Recommended Dilutions	Western Blot, Flow Cytometry, ELISA, Immunocytochemistry/ Immunofluorescence
Application Notes	Optimal dilution of this antibody should be experimentally determined. For optimal results using our Tandem dyes, please avoid prolonged exposure to light or extreme temperature fluctuations. These can lead to irreversible degradation or decoupling. When staining intracellular targets, specific attention to the fixation and permeabilization steps in your flow protocol may be required. Please contact our technical support team at technical@novusbio.com if you have any questions.

Images

Product Image: PD-1 Antibody (4F12) [Allophycocyanin/Cy7] [NBP1-04328APCCY7] - Vial of APC/Cy7 conjugated antibody. APC/Cy7 is optimally excited at 652 nm by the Red laser (633 or 640 nm) and has an emission maximum of 775 nm.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-04328APCCY7

NBP1-43317APCCY7	Mouse IgG2b Kappa Light Chain Isotype Control (MG2b) [Allophycocyanin/Cy7]
NBP1-77276PEP	PD-1 Antibody Blocking Peptide
210-TA-005	TNF-alpha [Unconjugated]
1086-PD-050	PD-1 [Unconjugated]

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-04328APCCY7

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications



