

Product Datasheet

UbcH8/Ube2L6 Antibody (k1H3) - BSA Free NBP1-04351

Unit Size: 0.1 ml

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-04351

Updated 9/9/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-04351



NBP1-04351

UbcH8/Ube2L6 Antibody (k1H3) - BSA Free


Product Information	
Unit Size	0.1 ml
Concentration	1.0 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	k1H3
Preservative	0.02% Sodium Azide
Isotype	IgG2b Kappa
Purity	Protein G purified
Buffer	PBS (pH 7.4), 10% Glycerol

Product Description	
Description	Novus Biologicals Mouse UbcH8/Ube2L6 Antibody (k1H3) - BSA Free (NBP1-04351) is a monoclonal antibody validated for use in IHC, WB, ELISA and ICC/IF. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Mouse
Gene ID	9246
Gene Symbol	UBE2L6
Species	Human
Immunogen	Recombinant human UbcH8/Ube2L6 (1-152aa) purified from E. coli


Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, ELISA, Immunocytochemistry/Immunofluorescence, Immunohistochemistry
Recommended Dilutions	Western Blot 1:2000, ELISA 1:100-1:2000, Immunohistochemistry 1:10-1:500, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry-Paraffin 1:10-1:500

Images


Western Blot: UbcH8/Ube2L6 Antibody (k1H3) [NBP1-04351] - The Cell lysate (20ug) were resolved by SDS-PAGE, transferred to PVDF membrane and probed with anti-human UBE2L6 antibody (1:2000). Proteins were visualized using a goat anti-mouse secondary antibody conjugated to HRP and an ECL detection system. Lane 1.: Jurkat cell lysate




Immunocytochemistry/Immunofluorescence: UbchH8/Ube2L6 Antibody (k1H3) [NBP1-04351] - ICC/IF analysis of UBE2L6 in HeLa cells line, stained with DAPI (Blue) for nucleus staining and monoclonal anti-human UBE2L6 antibody (1:100) with goat anti-mouse IgG-Alexa fluor 488 conjugate (Green).




Immunohistochemistry-Paraffin: UbchH8/Ube2L6 Antibody (k1H3) [NBP1-04351] - Ube2L6 Antibody (k1H3) [NBP1-04351]-Paraffin embedded sections of human breast lobule tissue were incubated with anti-human UBE2L6 for 2 hours at room temperature. Antigen retrieval was performed in 0.1M sodium citrate buffer and detected using Diaminobenzidine (DAB)



Western Blot: UbchH8/Ube2L6 Antibody (k1H3) [NBP1-04351] - Ube2L6 Antibody (k1H3) [NBP1-04351]-Cell lysates of MCF7 (20ug) were resolved by SDS-PAGE, transferred to NC membrane and probed with anti-human UBE2L6. Proteins were visualized using a goat anti-mouse secondary antibody conjugated to HRP and an ECL detection system.



Western Blot: UbchH8/Ube2L6 Antibody (k1H3) [NBP1-04351] - The Cell lysates (5ug) were resolved by SDS-PAGE, transferred to PVDF membrane and probed with anti-human UBE2L6 antibody (1:2000). Proteins were visualized using a goat anti-mouse secondary antibody conjugated to HRP and an ECL detection system. Lane 1.: 293T cell lysate Lane 2.: UBE2L6 Transfected 293T cell lysate





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-04351

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-43317-0.5mg	Mouse IgG2b Kappa Light Chain Isotype Control (MG2b)
NBP1-86885PEP	Ubch8/Ube2L6 Recombinant Protein Antigen

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-04351

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

