

Product Datasheet

G3BP1 Antibody NBP1-18922

Unit Size: 0.1 ml

Store at 4C. Do not freeze.

www.novusbio.com



technical@novusbio.com

Reviews: 1 Publications: 4

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-18922

Updated 9/9/2025 v.20.1

**Earn rewards for product
reviews and publications.**

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-18922



NBP1-18922**G3BP1 Antibody**

Product Information	
Unit Size	0.1 ml
Concentration	0.2 mg/ml
Storage	Store at 4C. Do not freeze.
Clonality	Polyclonal
Preservative	0.09% Sodium Azide
Isotype	IgG
Purity	Immunogen affinity purified
Buffer	TBS and 0.1% BSA


Product Description	
Description	Novus Biologicals Rabbit G3BP1 Antibody (NBP1-18922) is a polyclonal antibody validated for use in WB, ICC/IF and IP. Anti-G3BP1 Antibody: Cited in 3 publications. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	10146
Gene Symbol	G3BP1
Species	Human, Mouse, Mammal
Reactivity Notes	Mammal reactivity reported in scientific literature (PMID: 28794026).
Immunogen	The immunogen recognized by this antibody maps to a region between residue 275 and 325 of human GTPase activating protein (SH3 domain) binding protein 1 using the numbering given in entry NP_005745.1 (GeneID 10146).

Product Application Details	
Applications	Western Blot, Immunocytochemistry/ Immunofluorescence, Immunoprecipitation
Recommended Dilutions	Western Blot 1:2000-1:10000, Immunocytochemistry/ Immunofluorescence Validated from a verified customer review., Immunoprecipitation 2-10 ug/mg lysate




Images


Western Blot: G3BP1 Antibody [NBP1-18922] - Whole cell lysate (50 ug) from HeLa, MCF-7, U2OS, TCMK-1, and NIH 3T3 cells prepared using NETN lysis buffer. Antibody: Affinity purified rabbit anti-G3BP1 antibody used for WB at 0.04 ug/ml. Detection: Chemiluminescence with an exposure time of 1 second.



Immunocytochemistry/Immunofluorescence: G3BP1 Antibody [NBP1-18922] - G3BP1 in red (Alexa fluor 555) and nucleus in blue (DAPI). Human primary fibroblasts untreated in the left and treated in the right with G3BP1 stress granules labeling. Alexa fluor 555 1s exposure and DAPI 1s exposure. Fixation: PFA 4% 10min RT, permeabilization: Triton 0.1% TBS, blocking: Donkey Serum 10% TBS-T, Ab dilution: 1:1000, Secondary Ab: Alexa fluor 555 Donkey anti-rabbit 1:600, Mounting medium with DAPI. Image submitted by a verified customer review.



Western Blot: G3BP1 Antibody [NBP1-18922] - Whole cell lysate from HeLa (5, 15 and 50 mcg for WB; 1 mg for IP, 20% of IP loaded), 293T (T; 50 mcg) and mouse NIH3T3 (M; 50mcg) cells. Affinity purified rabbit anti-G3BP1 antibody used for WB at 0.04 mcg/ml (A) and 1 mcg/ml (B) and used for IP at 3 mcg/mg lysate. G3BP1 was also immunoprecipitated by rabbit anti-G3BP1 antibody NBP1-18923, which recognizes a downstream epitope.



Publications

Honda M, Oki S, Kimura R Et Al. High-depth spatial transcriptome analysis by photo-isolation chemistry Nature communications 2021-07-20 [PMID: 34285220] (IF/IHC)

Bussiere LD Examination of viral infection and host response to define therapeutic potential of oncolytic mammalian orthoreovirus Thesis 2020-01-01 (ICC/IF, Human)

Choudhury P, Bussiere L, Miller CL. Mammalian Orthoreovirus Factories Modulate Stress Granule Protein Localization by Interaction with G3BP1. J. Virol. 2017-08-09 [PMID: 28794026] (Mammal)

Fujimura K, Sasaki AT, Anderson P. Selenite targets eIF4E-binding protein-1 to inhibit translation initiation and induce the assembly of non-canonical stress granules. Nucleic Acids Res. 2012-06-20 [PMID: 22718973]



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-18922

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NBP1-50925-50ug	Recombinant Human G3BP1 His Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-18922

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

