

Product Datasheet

Cav3.2 Antibody (S55/10) - BSA Free NBP1-22444

Unit Size: 0.1 mg

Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Reviews: 1 Publications: 10

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-22444

Updated 2/21/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-22444



NBP1-22444

Cav3.2 Antibody (S55/10) - BSA Free

Product Information

Unit Size	0.1 mg
Concentration	1 mg/ml
Storage	Store at 4C short term. Aliquot and store at -20C long term. Avoid freeze-thaw cycles.
Clonality	Monoclonal
Clone	S55/10
Preservative	0.09% Sodium Azide
Isotype	IgG1
Purity	Protein G purified
Buffer	PBS (pH 7.4), 50% Glycerol

Product Description


Host	Mouse
Gene ID	8912
Gene Symbol	CACNA1H
Species	Human, Mouse, Rat
Specificity/Sensitivity	Detects approx 260 kDa. No cross-reactivity against Cav1.3.
Immunogen	Fusion protein amino acids 1019-1293 (II-III loop) of human Cav3.2

Product Application Details


Applications	Western Blot, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Frozen, Immunohistochemistry-Paraffin, Immunoprecipitation, Microarray
Recommended Dilutions	Western Blot 1:1000, Immunohistochemistry 1:1000, Immunocytochemistry/ Immunofluorescence 1:100, Immunoprecipitation, Immunohistochemistry-Paraffin 1:10-1:500, Immunohistochemistry-Frozen 1:10-1:500, Microarray
Application Notes	1 ug/ml of Cav3.2 Antibody was sufficient for detection of Cav3.2 in 10 ug of HEK cell lysate expressing Cav3.2 by colorimetric immunoblot analysis using Goat anti-mouse IgG:HRP as the secondary Antibody.

Images


Western Blot: Cav3.2 Antibody (S55/10) [NBP1-22444] - Western Blot analysis of Rat brain membrane lysate (native) showing detection of ~261 kDa Cav3.2 protein using Mouse Anti-Cav3.2 Monoclonal Antibody, Clone N55/10 (NBP1-22444). Block: 2% Skim Milk + 2% BSA in TBST. Primary Antibody: Mouse Anti-Cav3.2 Monoclonal Antibody (NBP1-22444) at 1:1000 for 2 hours at RT. Secondary Antibody: Anti-Mouse: HRP at 1:4000. Predicted/Observed Size: ~261 kDa.




Immunohistochemistry-Paraffin: Cav3.2 Antibody (S55/10) [NBP1-22444] - Tissue: hippocampus. Species: Human. Fixation: Bouin's Fixative and paraffin-embedded. Primary Antibody: Mouse Anti-CaV3.2 Calcium Channel Monoclonal Antibody at 1:1000 for 1 hour at RT. Secondary Antibody: FITC Goat Anti-Mouse (green) at 1:50 for 1 hour at RT.




Western Blot: Cav3.2 Antibody (S55/10) [NBP1-22444] - Western Blot using rat brain membranes and NBP1-22444 to target Cav3.2 channels visible at about 260 kDa.




Western Blot: Cav3.2 Antibody (S55/10) [NBP1-22444] - Western blot analysis of CaV3.2 proteins in mesenteric arteries of young versus old mice. (h), quantification of Western blot results. Mesenteric arteries were taken from 9 mice in each group. *, $p < .05$. Image collected and cropped by CiteAb from the following publication (onlinelibrary.wiley.com/doi/abs/10.1111/ace.13134) licensed under a CC-BY license.



Immunohistochemistry-Frozen: Cav3.2 Antibody (S55/10) [NBP1-22444] - Tissue: frozen brain section. Species: human. Fixation: 10% Formalin Solution for 12-24 hours at RT. Primary Antibody: Mouse Anti-CaV3.2 Calcium channel Monoclonal Antibody at 1:1000 for 1 hour at RT. Secondary Antibody: HRP/DAB Detection System: Biotinylated Goat Anti-Mouse, Streptavidin Peroxidase, DAB Chromogen (brown) for 30 minutes at RT. Counterstain: Mayer Hematoxylin (purple/blue) nuclear stain at 250-500 μ L for 5 minutes at RT.



Western Blot: Cav3.2 Antibody (S55/10) [NBP1-22444] - Age attenuates the role of CaV3.2 channels in Ca²⁺ spark generation & decreases CaV3.2 protein expression in VSMC. (a), Ca²⁺ fluorescence images of a Fluo-4-AM-loaded VSMC from a young mouse & time course of Ca²⁺ fluorescence changes in the cellular ROI (upper panel). Cell boundary is marked with dashed line. (b), same as (a) but in the presence of Ni²⁺ (50 μ M). (c), same as (a) but in a VSMC from an old mouse. (d), same as (c) but in the presence of Ni²⁺ (50 μ M). (e, f), summary of the results. Ca²⁺ spark frequency (e) & fraction of cells producing Ca²⁺ sparks (f) in VSMCs from young mice (n = 102), in VSMCs from young mice cells incubated with Ni²⁺ (n = 85), in VSMCs from aged mice (n = 129), & in VSMCs from aged mice cells incubated with Ni²⁺ (n = 127). Cells were isolated from 4 mice in each group; 25–40 cells were recorded & analyzed from each mouse. VSMC, vascular smooth muscle cell. (g), Western blot analysis of CaV3.2 proteins in mesenteric arteries of young versus old mice. (h), quantification of Western blot results. Mesenteric arteries were taken from 9 mice in each group. *, p < .05. n.s., not significant Image collected & cropped by CiteAb from the following publication (<https://pubmed.ncbi.nlm.nih.gov/32187825>), licensed under a CC-BY license. Not internally tested by Novus Biologicals.



Publications

- Cai H, Chen S, Sun Y et al. Interleukin-22 receptor 1-mediated stimulation of T-type Ca²⁺ channels enhances sensory neuronal excitability through the tyrosine-protein kinase Lyn-dependent PKA pathway. Cell communication and signaling : CCS 2024-06-04 [PMID: 38831315]
- Gandini MA, Souza IA, Khullar A Et al. Regulation of CaV3.2 channels by the receptor for activated C kinase 1 (Rack-1) Pflugers Archiv : European journal of physiology 2021-10-08 [PMID: 34623515] (WB)
- Zhang Y, Qian Z, Jiang D Et al. Neuromedin B receptor stimulation of Cav3.2 T-type Ca²⁺ channels in primary sensory neurons mediates peripheral pain hypersensitivity Theranostics 2021-09-09 [PMID: 34646374] (WB, Mouse)
- Fan G, Kamann M, Cui Y et Al. Age attenuates the T-type CaV 3.2-RyR axis in vascular smooth muscle Aging Cell 2020-03-18 [PMID: 32187825] (WB, Mouse)
- Wang H, Wei Y, Pu Y et al. Brain-derived neurotrophic factor stimulation of T-type Ca²⁺ channels in sensory neurons contributes to increased peripheral pain sensitivity Sci Signal 2019-09-24 [PMID: 31551295] (WB, Rat)
- Hashad AM, Harraz OF, Brett SE. Caveolae Link CaV3.2 Channels to BKCa-Mediated Feedback in Vascular Smooth Muscle Vascular Biology [PMID: 30354206] (WB)
- Garcia-Caballero A, Gandini MA, Huang S et al. Cav3.2 calcium channel interactions with the epithelial sodium channel ENaC Mol Brain 2019-02-08 [PMID: 30736831] (WB, Mouse)
- Garcia-Caballero A, Zhang FX, Hodgkinson V et al. T-type calcium channels functionally interact with spectrin (a/B) and ankyrin B Mol Brain 2018-05-02 [PMID: 29720258] (ICC/IF)
- Stemkowski P, Garcia-Caballero A, Gadotti VM et al. TRPV1 Nociceptor Activity Initiates USP5/T-type Channel-Mediated Plasticity. Cell Rep. 2016-12-13 [PMID: 27974205] (Mouse)
- Garcia-Caballero A, Gadotti VM, Stemkowski P et al. The Deubiquitinating Enzyme USP5 Modulates Neuropathic and Inflammatory Pain by Enhancing Cav3.2 Channel Activity. Neuron. 2014-09-03 [PMID: 25189210] (WB, IP, Mouse)



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-22444

HAF007	Goat anti-Mouse IgG Secondary Antibody [HRP]
NB720-B	Rabbit anti-Mouse IgG (H+L) Secondary Antibody [Biotin]
NBP1-97005-0.5mg	Mouse IgG1 Isotype Control (MG1)
NBP1-88176PEP	Cav3.2 Recombinant Protein Antigen

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-22444

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

