

# Product Datasheet

## CD38 Antibody (NIMR-5) [PE]

### NBP1-27954

Unit Size: 0.1 mg

Store at 4C. Do not freeze.

[www.novusbio.com](http://www.novusbio.com)



[technical@novusbio.com](mailto:technical@novusbio.com)

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:  
[www.novusbio.com/NBP1-27954](http://www.novusbio.com/NBP1-27954)

Updated 10/23/2024 v.20.1

Earn rewards for product  
reviews and publications.

Submit a publication at [www.novusbio.com/publications](http://www.novusbio.com/publications)

Submit a review at [www.novusbio.com/reviews/destination/NBP1-27954](http://www.novusbio.com/reviews/destination/NBP1-27954)



**NBP1-27954**

CD38 Antibody (NIMR-5) [PE]

Product Information	
Unit Size	0.1 mg
Concentration	0.1 mg/ml
Storage	Store at 4C. Do not freeze.
Clonality	Monoclonal
Clone	NIMR-5
Preservative	Sodium Azide
Isotype	IgG2a Kappa
Conjugate	PE
Purity	Protein A or G purified
Buffer	PBS and a stabilizing agent


Product Description	
Host	Rat
Gene ID	952
Gene Symbol	CD38
Species	Mouse
Specificity/Sensitivity	Mouse CD38
Immunogen	This CD38 Antibody (NIMR-5) was developed against BCL1 plasma membrane glycoproteins.
Notes	Protect conjugated forms from light.

Product Application Details	
Applications	Flow Cytometry
Recommended Dilutions	Flow Cytometry &lt;= 0.2 ug / 10 <sup>6</sup> cells




## Images

Flow Cytometry: CD38 Antibody (NIMR-5) [PE] [NBP1-27954] - BALB/c mouse splenocytes were stained with Rat Anti-Mouse CD38-PE and Rat Anti-Mouse CD45R-FITC.



Flow Cytometry: CD38 Antibody (NIMR-5) [PE] [NBP1-27954] - Analysis using the PE conjugate of NBP1-27951. Double staining of 0.3 ug/1<sup>6</sup> cells from BALB/c spleen with rat anti-mouse CD38-R-PE and rat anti-mouse CD19-FITC (Clone 6D5).





### **Novus Biologicals USA**

10730 E. Briarwood Avenue  
Centennial, CO 80112  
USA  
Phone: 303.730.1950  
Toll Free: 1.888.506.6887  
Fax: 303.730.1966  
nb-customerservice@bio-techne.com

### **Bio-Techne Canada**

21 Canmotor Ave  
Toronto, ON M8Z 4E6  
Canada  
Phone: 905.827.6400  
Toll Free: 855.668.8722  
Fax: 905.827.6402  
canada.inquires@bio-techne.com

### **Bio-Techne Ltd**

19 Barton Lane  
Abingdon Science Park  
Abingdon, OX14 3NB, United Kingdom  
Phone: (44) (0) 1235 529449  
Free Phone: 0800 37 34 15  
Fax: (44) (0) 1235 533420  
info.EMEA@bio-techne.com

### **General Contact Information**

www.novusbio.com  
Technical Support: nb-technical@bio-techne.com  
Orders: nb-customerservice@bio-techne.com  
General: novus@novusbio.com

### **Products Related to NBP1-27954**

---

NBP1-43780-0.1mg	Rat IgG2a Kappa Light Chain Isotype Control (R2a) [PE]
NBP1-27952	CD38 Antibody (NIMR-5) [FITC]
NBP2-52339-0.05mg	Recombinant Human CD38 His Protein
210-TA-005	TNF-alpha [Unconjugated]

---

### **Limitations**

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit [www.novusbio.com/guarantee](http://www.novusbio.com/guarantee)

Earn gift cards/discounts by submitting a review: [www.novusbio.com/reviews/submit/NBP1-27954](http://www.novusbio.com/reviews/submit/NBP1-27954)

Earn gift cards/discounts by submitting a publication using this product:  
[www.novusbio.com/publications](http://www.novusbio.com/publications)

