

Product Datasheet

CD3 epsilon Antibody (PPT3) [mFluor Violet 610 SE] NBP1-28225MFV610

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-28225MFV610

Updated 11/8/2023 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-28225MFV610



NBP1-28225MFV610

CD3 epsilon Antibody (PPT3) [mFluor Violet 610 SE]

| Product Information | |
|---------------------|---|
| Unit Size | 0.1 ml |
| Concentration | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage | Store at 4C in the dark. |
| Clonality | Monoclonal |
| Clone | PPT3 |
| Preservative | 0.05% Sodium Azide |
| Isotype | IgG1 Kappa |
| Conjugate | mFluor Violet 610 SE |
| Purity | Protein A or G purified |
| Buffer | 50mM Sodium Borate |

| Product Description | |
|-------------------------|--|
| Host | Mouse |
| Gene ID | 916 |
| Gene Symbol | CD3E |
| Species | Porcine |
| Specificity/Sensitivity | Porcine CD3 epsilon, Mr 20 kDa |
| Immunogen | Purified CD3 molecules from porcine thymus |

| Product Application Details | |
|-----------------------------|--|
| Applications | Flow Cytometry |
| Recommended Dilutions | Flow Cytometry |
| Application Notes | Optimal dilution of this antibody should be experimentally determined. |

Images

CD3 epsilon Antibody (PPT3) [mFluor Violet 610 SE] [NBP1-28225MFV610] - Vial of mFluor Violet 610 conjugated antibody. mFluor Violet 610 is optimally excited at 421 nm by the Violet laser (405 nm) and has an emission maximum of 613 nm.



mFluor™ Violet 610

| LASER (nm) | FILTER |
|---------------------|-------------------|
| Violet (405) | 605/30 |
| EXCITATION MAX (nm) | EMISSION MAX (nm) |
| 421 | 613 |



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-28225MFV610

| | |
|-------------|---|
| NBP2-22752 | Recombinant Human CD3 epsilon His Protein |
| 210-TA-005 | TNF-alpha [Unconjugated] |
| 9850-CD-050 | CD3 epsilon [Unconjugated] |
| 285-IF-100 | IFN-gamma [Unconjugated] |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-28225MFV610

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

