

Product Datasheet

Importin alpha 3/KPNA4 Antibody NBP1-31260

Unit Size: 0.1 ml

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Reviews: 1 Publications: 4

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-31260

Updated 10/23/2024 v.20.1

**Earn rewards for product
reviews and publications.**

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-31260



NBP1-31260

Importin alpha 3/KPNA4 Antibody


Product Information	
Unit Size	0.1 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.01% Thimerosal
Isotype	IgG
Purity	Antigen Affinity-purified
Buffer	PBS, 1% BSA, 20% Glycerol (pH7)
Target Molecular Weight	58 kDa

Product Description	
Host	Rabbit
Gene ID	3840
Gene Symbol	KPNA4
Species	Human, Mouse, Rat
Immunogen	Carrier-protein conjugated synthetic peptide encompassing a sequence within the C-terminus region of human Importin alpha 3/KPNA4. The exact sequence is proprietary.


Product Application Details	
Applications	Western Blot, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin
Recommended Dilutions	Western Blot 1:500-1:3000, Immunohistochemistry 1:100-1:1000, Immunocytochemistry/ Immunofluorescence 1:100-1:1000, Immunohistochemistry-Paraffin 1:100-1:1000

Images


Western Blot: Importin alpha 3/KPNA4 Antibody [NBP1-31260] - Sample (30 ug of whole cell lysate) A: HeLa B: Hep G2 7.5% SDS PAGE diluted at 1:1000




Immunocytochemistry/Immunofluorescence: Importin alpha 3/KPNA4 Antibody [NBP1-31260] - Analysis of HeLa, using KPNA4 antibody at 1:200 dilution.




Immunohistochemistry-Paraffin: Importin alpha 3/KPNA4 Antibody [NBP1-31260] - H1299 xenograft, using KPNA4 antibody at 1:500 dilution. Antigen Retrieval: Trilogy™ (EDTA based, pH 8.0) buffer, 15min.




Western Blot: Importin alpha 3/KPNA4 Antibody [NBP1-31260] - Sample (30 ug of whole cell lysate) A: 293T B: HepG2 C: Molt-4 7.5% SDS PAGE diluted at 1:150



Western Blot: Importin alpha 3/KPNA4 Antibody [NBP1-31260] - A. 30 ug Rat2 whole cell lysate/extract 7.5 % SDS-PAGE KPNA4 antibody [C3], C-term dilution: 1:1000



Western Blot: Importin alpha 3/KPNA4 Antibody [NBP1-31260] - A. 30 ug Neuro2A whole cell lysate/extract. B. 30 ug GL261 whole cell lysate/extract. C. 30 ug C8D30 whole cell lysate/extract. D. 30 ug NIH-3T3 whole cell lysate/extract. E. 30 ug BCL-1 whole cell lysate/extract. 7.5 % SDS-PAGE KPNA4 antibody [C3], C-term dilution: 1:1000



Publications

Andreu I, Granero-Moya I, Chahare NR et al. Mechanical force application to the nucleus regulates nucleocytoplasmic transport *Nature cell biology* 2022-06-01 [PMID: 35681009] (Mouse)

Solomon DA, Stepto A, Au WH et al. A feedback loop between dipeptide-repeat protein, TDP-43 and karyopherin-? mediates C9orf72-related neurodegeneration. *Brain* 2018-10-01 [PMID: 30239641] (WB, IF/IHC, ICC/IF, Human, Mouse)

Yang J, Lu C, Wei J et al. Inhibition of KPNA4 attenuates prostate cancer metastasis. *Oncogene* 2017-01-01 [PMID: 27941876] (WB, Mouse)

Jeong SA, Kim K, Lee JH et al. Akt-mediated phosphorylation increases the binding affinity of hTERT for importin alpha to promote nuclear translocation. *J. Cell. Sci.* 2015-05-21 [PMID: 25999477] (WB, Human)



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-31260

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
H00003840-P01-10ug	Recombinant Human Importin alpha 3/KPNA4 GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-31260

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

