

Product Datasheet

Enolase 1 Antibody - BSA Free NBP1-31324

Unit Size: 100 ul

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-31324

Updated 2/21/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-31324



NBP1-31324

Enolase 1 Antibody - BSA Free

Product Information

Unit Size	100 ul
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.01% Thimerosal
Isotype	IgG
Purity	Antigen Affinity-purified
Buffer	0.1M Tris (pH 7), 0.1M Glycine, 10% Glycerol
Target Molecular Weight	47 kDa

Product Description


Host	Rabbit
Gene ID	2023
Gene Symbol	ENO1
Species	Human, Mouse, Rat, Zebrafish
Reactivity Notes	Xenopus laevis (90%).
Immunogen	Recombinant protein encompassing a sequence within the center region of human Enolase 1. The exact sequence is proprietary.

Product Application Details


Applications	Western Blot, Dot Blot, Immunocytochemistry/ Immunofluorescence, Immunohistochemistry, Immunohistochemistry-Paraffin
Recommended Dilutions	Western Blot 1:500-1:20000, Immunohistochemistry 1:100-1:1000, Immunocytochemistry/ Immunofluorescence 1:100-1:1000, Immunohistochemistry-Paraffin 1:100-1:1000, Dot Blot
Application Notes	Dot Blot-Assay dependent

Images


Western Blot: Enolase 1 Antibody [NBP1-31324] - Various whole cell extracts (30 ug) were separated by 10% SDS-PAGE, and the membrane was blotted with ENO1 antibody [N3C3] diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody (NBP2-19301) was used to detect the primary antibody.




Immunocytochemistry/Immunofluorescence: Enolase 1 Antibody [NBP1-31324] - HeLa cells were fixed in 4% paraformaldehyde at RT for 15 min. Green: ENO1 protein stained by ENO1 antibody [N3C3] diluted at 1:200. Red: alpha Tubulin, a cytoskeleton marker, stained by alpha Tubulin antibody [B-5-1-2] diluted at 1:10000. Blue: Hoechst 33342 staining.




Immunohistochemistry-Paraffin: Enolase 1 Antibody [NBP1-31324] - Paraffin-embedded gastric tissue, using antibody at 1:100 dilution.




Western Blot: Enolase 1 Antibody [NBP1-31324] - Sample (50 ug of whole cell lysate) A: Mouse brain 10% SDS PAGE; antibody diluted at 1:1000.




Western Blot: Enolase 1 Antibody [NBP1-31324] - Sample (30 ug of whole cell lysate) A: A431 B: H1299 10% SDS PAGE, antibody diluted at 1:3000.



Western Blot: Enolase 1 Antibody [NBP1-31324] - Non-transfected (-) and transfected (+) 293T whole cell extracts (30 ug) were separated by 10% SDS-PAGE, and the membrane was blotted with ENO1 antibody [N3C3].



Western Blot: Enolase 1 Antibody [NBP1-31324] - Rat tissue extract (50 ug) was separated by 10% SDS-PAGE, and the membrane was blotted with ENO1 antibody [N3C3] diluted at 1:1000. The HRP-conjugated anti-rabbit IgG antibody (NBP2-19301) was used to detect the primary antibody, and the signal was developed with Trident ECL plus-Enhanced.




Western Blot: Enolase 1 Antibody [NBP1-31324] - Sample (30 ug of whole cell lysate)


A: zebrafish eye

10% SDS PAGE

NBP1-31324 diluted at 1:1000




Immunohistochemistry: Enolase 1 Antibody [NBP1-31324] - Immunohistochemical analysis of paraffin-embedded zebrafish tissue, using Enolase 1 antibody [N3C3] (NBP1-31324) at 1:300 dilution.



Immunohistochemistry: Enolase 1 Antibody [NBP1-31324] - Enolase 1 antibody [N3C3] detects Eno1 protein on zebrafish by whole mount immunohistochemical analysis.

Sample: 2 days-post-fertilization zebrafish embryo.

Enolase 1 antibody [N3C3] (NBP1-31324) dilution: 1:100.





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-31324

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
NBP1-30262	Recombinant Human Enolase 1 Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-31324

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

