

Product Datasheet

PPAR alpha/NR1C1 Antibody - BSA Free NBP1-31834

Unit Size: 0.1 ml

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-31834

Updated 2/21/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-31834



NBP1-31834

PPAR alpha/NR1C1 Antibody - BSA Free

Product Information

Unit Size	0.1 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.025% Proclin 300
Isotype	IgG
Purity	Antigen Affinity-purified
Buffer	PBS, 20% Glycerol
Target Molecular Weight	52 kDa

Product Description


Host	Rabbit
Gene ID	5465
Gene Symbol	PPARA
Species	Human, Mouse, Rat, Insect, Plant
Reactivity Notes	Queen Bee, Honeybee. Bovine (88%), Guinea pig (87%)
Immunogen	Recombinant protein encompassing a sequence within the N-terminus region of human PPAR alpha/NR1C1. The exact sequence is proprietary.

Product Application Details


Applications	Western Blot, Immunohistochemistry, Immunohistochemistry-Paraffin, Chromatin Immunoprecipitation (ChIP)
Recommended Dilutions	Western Blot 1:500-1:3000, Immunohistochemistry, Immunohistochemistry-Paraffin Assay dependent, Chromatin Immunoprecipitation (ChIP) Assay dependent

Images


Western Blot: PPAR alpha/NR1C1 Antibody [NBP1-31834] - Various tissue extracts (50 ug) were separated by 10% SDS-PAGE, and the membrane was blotted with PPAR alpha antibody diluted at 1:500. The HRP-conjugated anti-rabbit IgG antibody (NBP2-19301) was used to detect the primary antibody.




Chromatin Immunoprecipitation: PPAR alpha/NR1C1 Antibody [NBP1-31834] - Cross-linked ChIP was performed with HepG2 chromatin extract and 5 ug of either control rabbit IgG or anti-PPAR alpha antibody. The precipitated DNA was detected by PCR with primer set targeting to C3 promotor or NOD1 gene.




Immunohistochemistry-Paraffin: PPAR alpha/NR1C1 Antibody [NBP1-31834] - PPAR alpha/NR1C1 antibody detects PPAR alpha/NR1C1 protein at nucleus by immunohistochemical analysis. Sample: Paraffin-embedded rat brain. PPAR alpha/NR1C1 stained by PPAR alpha/NR1C1 antibody (NBP1-31834) diluted at 1:165. Antigen Retrieval: Citrate buffer, pH 6.0, 15 min



Western Blot: PPAR alpha/NR1C1 Antibody [NBP1-31834] - Whole cell extract (30 ug) was separated by 10% SDS-PAGE, and the membrane was blotted with PPAR alpha/NR1C1 antibody (NBP1-31834) diluted at 1:500. The HRP-conjugated anti-rabbit IgG antibody was used to detect the primary antibody.



Immunohistochemistry-Paraffin: PPAR alpha/NR1C1 Antibody [NBP1-31834] - PPAR alpha/NR1C1 antibody detects PPAR alpha/NR1C1 protein at nucleus by immunohistochemical analysis. Sample: Paraffin-embedded mouse intestine. PPAR alpha/NR1C1 stained by PPAR alpha/NR1C1 antibody (NBP1-31834) diluted at 1:165. Antigen Retrieval: Citrate buffer, pH 6.0, 15 min





Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-31834

HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control
H00005465-P01-10ug	Recombinant Human PPAR alpha/NR1C1 GST (N-Term) Protein

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-31834

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

