

Product Datasheet

ORP150/HSP12A Antibody - BSA Free NBP1-32140

Unit Size: 0.1 ml

Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.

www.novusbio.com



technical@novusbio.com

Publications: 1

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NBP1-32140

Updated 9/25/2025 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NBP1-32140



NBP1-32140

ORP150/HSP12A Antibody - BSA Free


Product Information	
Unit Size	0.1 ml
Concentration	Concentrations vary lot to lot. See vial label for concentration. If unlisted please contact technical services.
Storage	Aliquot and store at -20C or -80C. Avoid freeze-thaw cycles.
Clonality	Polyclonal
Preservative	0.01% Thimerosal
Isotype	IgG
Purity	Antigen Affinity-purified
Buffer	0.1M Tris, 0.1M Glycine, 10% Glycerol
Target Molecular Weight	111 kDa

Product Description	
Description	Novus Biologicals Rabbit ORP150/HSP12A Antibody - BSA Free (NBP1-32140) is a polyclonal antibody validated for use in IHC, WB and ICC/IF. Anti-ORP150/HSP12A Antibody: Cited in 1 publication. All Novus Biologicals antibodies are covered by our 100% guarantee.
Host	Rabbit
Gene ID	10525
Gene Symbol	HYOU1
Species	Human, Mouse, Rat, Plant
Reactivity Notes	Bovine (88%).
Immunogen	Recombinant protein encompassing a sequence within the C-terminus region of human ORP150/HSP12A. The exact sequence is proprietary.


Product Application Details	
Applications	Western Blot, Immunohistochemistry-Paraffin, Immunocytochemistry/Immunofluorescence, Immunohistochemistry
Recommended Dilutions	Western Blot 1:500-1:3000, Immunohistochemistry 1:100-1:1000, Immunocytochemistry/ Immunofluorescence 1:100-1:1000, Immunohistochemistry-Paraffin 1:100-1:1000

Images


Western Blot: ORP150/HSP12A Antibody [NBP1-32140] - Mouse tissue extracts (50 ug) was separated by 5% SDS-PAGE, and blotted with ORP150 antibody diluted by 1:500. The HRP-conjugated anti-rabbit IgG antibody was used to detect the primary antibody.




Immunocytochemistry/Immunofluorescence: ORP150/HSP12A Antibody [NBP1-32140] - HeLa cells were fixed in 4% paraformaldehyde at RT for 15 min. Green: ORP150 protein stained by ORP150 antibody [C2C3], C-term diluted at 1:200. Red: phalloidin, a cytoskeleton marker, stained by phalloidin (invitrogen, A12380) diluted at 1:200. Blue: Hoechst 33342 staining.




Immunohistochemistry-Paraffin: ORP150/HSP12A Antibody [NBP1-32140] - Meta-analysis of significant metastatic upregulated proteins in acid exosomes through Prognoscan. Confocal immunofluorescence analysis with the indicated antibodies on primary (apical and mid region), and metastatic melanoma biopsies from the same patient. Ki67 labeling was used as internal control. Red, CFL, GSN, HYOU1; blue, DAPI; green, Ki67. Bar 50 um Image collected and cropped by CiteAb from the following publication (<https://jeccr.biomedcentral.com/articles/10.1186/s13046-018-0915-z>), licensed under a CC-BY license.




Western Blot: ORP150/HSP12A Antibody [NBP1-32140] - Various whole cell extracts (30 ug) were separated by 5% SDS-PAGE, and blotted with ORP150 antibody diluted by 1:1000. The HRP-conjugated anti-rabbit IgG antibody was used to detect the primary antibody.




Western Blot: ORP150/HSP12A Antibody [NBP1-32140] - Rat tissue extracts (50 ug) was separated by 5% SDS-PAGE, and blotted with ORP150 antibody diluted by 1:1000. The HRP-conjugated anti-rabbit IgG antibody was used to detect the primary antibody.




Immunocytochemistry/Immunofluorescence: ORP150/HSP12A Antibody [NBP1-32140] - Analysis of methanol-fixed Hep3B, using antibody at 1:500 dilution.




Costained with Hoechst 33342




Immunohistochemistry-Paraffin: ORP150/HSP12A Antibody [NBP1-32140] - Paraffin-embedded CL1-0 xenograft, using antibody at 1:100 dilution.



Immunohistochemistry-Paraffin: ORP150/HSP12A Antibody [NBP1-32140] - Paraffin-embedded mouse fore brain. ORP150 antibody [C2C3], C-term diluted at 1:500.



Immunocytochemistry/ Immunofluorescence: ORP150/HSP12A Antibody [NBP1-32140] - Meta-analysis of significant metastatic upregulated proteins in acid exosomes through Prognoscan. a Kaplan-Meyer curves depicting the Overall Survival of metastatic melanoma patients with high (red lines) or low (green lines) expression of the genes listed in Table 1 that resulted statistically significant. Mentel-Cox p-values for each plot are showed. b Western blotting analysis of ctr & pH 6.0 exosomes (20 μ g). Cell lysate was used as positive control. Expression of CFL, GSN & HSP90 $\alpha\beta$ is shown. Ponceau red is shown for equal exosome protein amount. c confocal immunofluorescence analysis with the indicated antibodies on primary (apical & mid region), & metastatic melanoma biopsies from the same patient. Ki67 labeling was used as internal control. Red, CFL, GSN, HYOU1; blue, DAPI; green, Ki67. Bar 50 μ m Image collected & cropped by CiteAb from the following publication (<https://pubmed.ncbi.nlm.nih.gov/30290833>), licensed under a CC-BY license. Not internally tested by Novus Biologicals.




ORP150/HSP12A antibody [C2C3], C-term detects ORP150/HSP12A protein at endoplasmic reticulum by immunofluorescent analysis.
Sample: HeLa cells were fixed in 4% paraformaldehyde at RT for 15 min.

Green: ORP150/HSP12A protein stained by ORP150/HSP12A antibody [C2C3], C-term (NBP1-32140) diluted at 1:200.

Red: phalloidin, a cytoskeleton marker, stained by phalloidin (invitrogen, A12380) diluted at 1:200.

Blue: Hoechst 33342 staining.




ORP150/HSP12A antibody [C2C3], C-term detects ORP150/HSP12A protein at cytoplasm by immunofluorescent analysis.

Sample: U2OS cells were fixed in 4% paraformaldehyde at RT for 15 min.

Green: ORP150/HSP12A protein stained by ORP150/HSP12A antibody [C2C3], C-term (NBP1-32140) diluted at 1:1000.

Red: alpha Tubulin, a cytoskeleton marker, stained by alpha Tubulin antibody [GT114] diluted at 1:1000.

Blue: Hoechst 33342 staining.



Publications

Boussadia Z, Lamberti J, et al. Acidic microenvironment plays a key role in human melanoma progression through a sustained exosome mediated transfer of clinically relevant metastatic molecules. J Exp Clin Cancer Res 2018-10-05 [PMID: 30290833] (IF/IHC, Human)



Novus Biologicals USA

10730 E. Briarwood Avenue
Centennial, CO 80112
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
nb-customerservice@bio-techne.com

Bio-Techne Canada

21 Canmotor Ave
Toronto, ON M8Z 4E6
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada.inquires@bio-techne.com

Bio-Techne Ltd

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info.EMEA@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: nb-technical@bio-techne.com
Orders: nb-customerservice@bio-techne.com
General: novus@novusbio.com

Products Related to NBP1-32140

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
HAF008	Goat anti-Rabbit IgG Secondary Antibody [HRP]
NB7160	Goat anti-Rabbit IgG (H+L) Secondary Antibody [HRP]
NBP2-24891	Rabbit IgG Isotype Control

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NBP1-32140

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

